IC10/20FRLCD

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IC10/20FRLCD

Auto-Reader Incubator

For incubation and readout of Rapid, Super Rapid or Ultra Rapid Biological Indicators and Hygiene Monitoring Systems

Usage

For incubation and readout of:

Bionova® Rapid Fluorescence Readout Biological Indicators: BT220 (Steam): 3 hs. at 60 °C, BT221 (Steam): 1 h. at 60 °C, BT110 (EO): 4 hs. at 37 °C, BT102 (Formaldehyde): 2 hs. at 60°C, BT95 (VH₂O₂): 2 hs. at 60°C.

Bionova® Super Rapid Fluorescence Readout Biological Indicators: BT222 (Steam): 1 h. at 60 °C, BT223 (Steam): 30 min. at 60 °C, BT96 (VH₂O₂): 30 min. at 60°C.

Bionova® Ultra Rapid Fluorescence Readout Biological Indicators: BT224 (Steam): 20 min. at 60°C.

Chemdye® PRO1 MICRO Protein Pen System: 15 min at 60 °C

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

Characteristics

12 metal positions (3.8 cm deep, 0.9 cm diameter) to incubate Rapid/Super Rapid/Ultra Rapid Fluorescence Readout Self-contained biological indicators

1 position (3.8 cm deep, 1.3 cm diameter) to incubate Protein Pen

1 ampoule crusher

1 special hole for external thermometer.

3.5 " touchscreen LCD. Visual information of the remaining time of the incubation process.

Temperature display on the screen and remaining time for each position.

Dual Temperature System: Allows selecting between two different incubation temperatures (37°C y 60°C).

Allows running different incubation times simultaneously, provided they are of equal incubation temperature.

Audible alarm for event indications.

ETHERNET connection for recording results on PC using a Reading and Traceability software.

Remote display of the screen in PC and Smartphones.

It stores the last 208 results.

Thermal Printer of Results. Allows reprinting the last 208 results.

Automatic detection and cancelation of biological indicator fluorescence reading.

Firmware update by the user.

Fluorescence reading wavelength: ~ 460 nm

Dimensions: 18.0 cm high, 26.0 cm diameter.

Voltage range: 100 – 240 V AC

Power: 28 W

Frequency: 50-60 Hz

Environmental conditions during manufacture

Room temperature 15-35 °C, RH 35-85 %

Storage conditions

T = 10-30 °C, RH 30-80 %, keep away from light.

Transportation conditions

Storage conditions should be strictly followed. The product should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf-life

5 years

Packing

1 unit per box.

Packaging information: product code and description, storage conditions and manufacturer information.

Labelling

Same label on product and product's box: product code and description, storage conditions, batch number, serial number, MAC number, bar code and manufacturer information.

Possible target markets

Healthcare, Food, Pharmaceutical and medical devices industries.

Other important information

The IC10/20FRLCD incubator has an icon on its screen that allows selecting the temperature at 37 °C or 60 °C.

In addition, it allows the configuration of each reading position with a different incubation program.

The IC10/20FRLCD incubator has an automatic readout system and a printer to record data.

It also has a history of records to query the last 208 results, which can be viewed from the device screen or remotely accessing the embedded web server.

Description of identified and/or predictable risks which could represent a risk for the product

Do not pour any liquid inside.

Do not immerse into any liquid.

Only for indoor use.

Disconnect the power cord before cleaning.

Do not use abrasive or corrosive cleaners or disinfectants.

Make sure that the incubator is connected to a properly rated power cord.

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