PCD220-2

Process Challenge Device

Test pack for Steam sterilization processes

Usage

For routine load release, especially for loads containing implants. For routine sterilizer monitoring. For periodic validation of the sterilizer (after installation, repair, relocation)

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity-air displacement steam autoclaves.

Conditions: Vacuum assisted steam sterilization processes and gravity-air displacement steam sterilization processes between 121 °C - 135 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11140-1:2014, ISO 11138-1:2017, ISO 11138-3:2017, ANSI/AAMI ST79

Classification

Class 1, according to the risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

FDA 510(k)

K163646

Characteristics

Porous cards system holding a Bionova® BT220 Self-contained Biological Indicator and a Type 5 Integrator Indicator printed on PCDBI-2-RC self-adhesive Record Card, where data on sterilization cycle may be written. Two possibilities for documentation record: the whole Record Card can be adhered, or alternatively, only the area containing the indicator may be recorded.

Cardboard box (10.5 cm x 6.5 cm x 1.8 cm approximately) with a printed Type 1 Process Indicator.

Process Indicator

Initial color: **blue**

Final color: **dark grey or black**

PCDBI-2-RC Integrator Indicator (within the challenge device)

Initial color: **purple**

Final color: **green**

100 % Toxic Heavy Metals free

Integration condition is calibrated against the kill time of a 10⁶ G. *stearothermophilus* ATCC 7953 spores' population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

BT220 Rapid Self-contained Biological Indicator (within the challenge device): $\geq 10^6$ *G.stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

Final results: by fluorescence after incubation for 3 hours at 60 °C (sensitivity $\geq 97\%$). An optional visual pH color change confirmation could be made after 48 hours of incubation. A 7-day readout is optional to confirm 3-hour results by color change.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

T =15-30 °C, RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place in its original box.

Transport conditions

Storage conditions should be strictly followed. Carry in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf life

2 years.

Packing

Packing information: product code and description, regulation, process for intended use, conditions and instructions for use, storage conditions, presentation, manufacturer data and label information.

Labelling

On the package: product code, process for intended use, batch number, batch number of the Self-Contained Biological Indicator, manufacture and expiration date, regulation, bar code and datamatrix.

On the BT220 Self-Contained Biological Indicator: label printed with Steam reactive ink (turns to brown), product code, batch number, expiration date, process for intended use and name of the organism.

On the Chemical Indicator Integrator PCDBI-2-RC: product code, process for intended use and classification according to the regulation.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

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