Bionova® Cloud Traceability Software

Automatic System for Quality Control &
Traceability of Infection Control Processes

Rev. 0 | April 2021

Download the latest version of the user´s manual and other resources for your equipment here.

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4 Access rights and switching users

4 User account management

4 Resetting and changing the password

4 Modify company settings

5 Company events registration

5 Language setting

5 Main screen

6 Presentation of the information in Bionova® Cloud

6 Selection or edition of multiple records

7 Storage of results

7 Areas

7 Users

7 Connected incubators

7 Biological indicators

7 In process

7 Finished

8 Biological indicators for sterilization - Manual data

8 Biological indicators for disinfection - Manual data

9 Calculation of the sterilization process risk (SPR)

11 Protein indicators

11 In progress

11 Finished

12 Information to be completed manually

12 HTM 01-01

13 Control analysis

13 Results interpretation

14 Chemical indicators

14 Detection

15 Information to be completed manually

15 Equipment

16 Licenses

16 My licenses

16 Filed results

16 Reports download

16 Equipment calibration

17 Other functions

17 Clock synchronization

17 Full-screen mode

■ Access


To access the Bionova® Cloud Traceability Software, open any web browser and enter the following URL: <https://bionovacloud.com>

Login

To access the system, you need your login credentials, an email, and a password.

You can create a new company to start operating with the system. Note that a new company has its own data set, it cannot be shared with an existing company.

If you do not remember your password, you can reset it from here by clicking on "Forgot your Password?"

■ Access rights and switching users

The Administrator user (Admin) has full access rights to the functioning of the application.

Menu item: COMPANY > USERS.

When using this user account, the Administrator can:

_Add or remove other user accounts (see *User Account Management*).

_Reset passwords for existing accounts (see *Resetting and changing the password*).

_Modify application settings (see *Modify Company settings*).

_Display the company's entire operating history (see *Company events registration*).

Additional user accounts (corresponding to the Operators) can only see the results of the analysis performed by the user of the account, but not those done by other users.

If different user accounts have been enabled, you can switch from one user to another by clicking on the user name in the "Logout" option.

Note: In case the Admin user forgets the password, click on *I forgot my password*, and an email will be sent to reset it.

User account management

_With the "Add" button, the Admin user can add a new user account to the system.

_By pointing the mouse cursor at the table where existing users are displayed, the Admin can edit user accounts, view user account details, disable existing accounts, and reset passwords.

Resetting and changing the password

The Admin user can reset the password of other users. This feature can be useful in case a user has forgotten his password. Once the administrator has reset a user's password, the user can enter the program using the self-generated password sent to him/her by email.

Also, all users can change their password using the option: CONFIGURATION > CHANGE PASSWORD.

Modify company settings

The Admin user can define the general settings of the program: CONFIGURATION > PARAMETERS.

RESET ALL DEFAULT VALUES	
Company Security Mode	
Relaxed	
<input type="button" value="SAVE"/>	<input type="button" value="RESET DEFAULT VALUE"/>
Company Language	
English	
<input type="button" value="SAVE"/>	<input type="button" value="RESET DEFAULT VALUE"/>
Days For License Expiration Notification	
30	
<input type="button" value="SAVE"/>	<input type="button" value="RESET DEFAULT VALUE"/>
Protein Umbra	
1	
<input type="button" value="SAVE"/>	<input type="button" value="RESET DEFAULT VALUE"/>

↳ Access rights and switching users

1 | The security mode defines the configuration that affects the amount of mandatory fields needed to file a new record, being "Relaxed" the configuration with the least amount of data to be completed.

2 | Change the company language.

3 | Set the number of days before the system notifies the expiration of a license.

4 | The Threshold defines the protein level (μg) at which the hygiene indicator result is considered positive. Default reference value: 1 μg .

This information is very important as it will be displayed on all reports. Once the configuration is complete, the Admin user can start operating the program.

Company events registration

The Admin user can monitor the movements that take place

within the company, using the option: CONFIGURATION > LOG. The user who registered a new movement, the exact date of its completion, and a detail of what he/she has done can be seen.

■ Language setting

The Admin user can modify the app language, as long as he/she takes into account that all users of that company will be reached by this change (see *Modify company settings*).



■ Main screen

The main screen of the app divides into three parts:

1 | **Header:** contains the current user's session information and a button to collapse the side menu.

2 | **Main menu:** a panel with graphics and statistics about Bionova® Cloud is displayed.

3 | **Side menu:** contains all available screens for navigation.

■ Presentation of the information in Bionova® Cloud

Bionova® Cloud presents a standardized way of representing information for all screens.

STERILIZER 1		
Q 2		
NAME ↑	BRAND	SERIAL
STERILIZER 1	STE 1	1111
STERILIZER 2	STE 2	2222
STERILIZER 3	STE 3	3333
STERILIZER 4	STE 4	4444

ADD

5

ROWS PER PAGE 50 ▾

(1-4/4) 6

1 | Title of the table.

2 | Search engine: it allows to filter the table records by

keyword. It filters through all the columns in the table.

3 | Header: it contains the name of each column, and allows you to sort the body records by clicking on it.

4 | Body: the table body contains the information on the records.

5 | Action buttons: the buttons that allow you to perform actions in the table.

6 | Table footer: it contains the pagination information.

a. The number of records shown per page.

b. The records being displayed on the current page, based on the total records contained in the table.

c. Browser buttons between pages.

7 | Actions panel for each record in the table. It is displayed only when the cursor is pointing to a specific register.

8 | When the table is wider than the screen, you can scroll horizontally by holding down the left click and moving the mouse to the sides to navigate to the next columns (the cursor will turn into a hand).

■ Selection or edition of multiple records

PROVISIONAL									
	READ NUMBER ↓	POSITION NUMBER	START TIME	END TIME	PROGRAM	RESULT	AVERAGE TEMPERATURE	TICKET NUMBER	INCUBATOR
<input checked="" type="checkbox"/>	209	9	12/17/2020 11:59:56	12/17/2020 12:19:58	ULTRA	POSITIVE	60,0 °C	204	IC1020FR_TERRA 2120-999
<input checked="" type="checkbox"/>	208	7	12/17/2020 11:59:56	12/17/2020 12:19:57	ULTRA	NEGATIVE	60,0 °C	203	IC1020FR_TERRA 2120-999
<input checked="" type="checkbox"/>	206	4	12/17/2020 11:53:42	12/17/2020 14:53:42	3 HS.	NEGATIVE	59,9 °C	206	IC1020FR_TERRA 2120-999
<input checked="" type="checkbox"/>	205	4	12/17/2020 11:53:39	12/17/2020 11:53:40	3 HS.	CANCELED	60,0 °C	202	IC1020FR_TERRA 2120-999
<input checked="" type="checkbox"/>	203	2	12/17/2020 11:18:39	12/17/2020 11:46:04	3 HS.	CANCELED	60,0 °C	201	IC1020FR_TERRA 2120-999
<input checked="" type="checkbox"/>	202	2	12/17/2020 11:18:22	12/17/2020 11:18:24	3 HS.	CANCELED	60,0 °C	199	IC1020FR_TERRA 2120-999
<input checked="" type="checkbox"/>	197	2	12/16/2020 11:34:54	12/16/2020 11:37:09	3 HS.	CANCELED	60,0 °C	196	IC1020FR_TERRA 2120-999
<input type="checkbox"/>	196	1	12/16/2020 11:34:56	12/16/2020 11:37:08	3 HS.	CANCELED	60,0 °C	195	IC1020FR_TERRA 2120-999
<input type="checkbox"/>	195	1	12/16/2020 11:31:31	12/16/2020 11:33:40	3 HS.	CANCELED	59,9 °C	194	IC1020FR_TERRA 2120-999
<input type="checkbox"/>	194	1	12/16/2020 10:47:54	12/16/2020 10:47:56	3 HS.	CANCELED	60,0 °C	192	IC1020FR_TERRA 2120-999

ROWS PER PAGE 10 ▾

(1-10/32)

> >|

2 MULTIPLE EDIT

MULTIPLE TO HISTORICAL 3

1 | Selection column: when you click on a record of the table, this column will appear with a checkbox to select multiple options.

2 | When some records are selected, the button for multiple editing will be available.

3 | It will also be available a button to "File" several records at a time.

■ Storage of results

When a set of results is stored, these will become part of "saved" readouts (see *Saved results*) being stored and available for future reference. Before the permanent storage of data, the program performs a series of verifications to ensure its integrity. It ensures that the user does not lose any relevant information about these tests when saving the result.

The mandatory data for result archiving depends on the "Company Security Mode" setting (see *Modify company settings*).

Storing an individual result

To permanently store the information corresponding to a particular indicator, open the action panel and click on to saved. (see *Action panel*).

Storing multiple results

To permanently store two or more results in one step (see *Multiple selections*).

■ Areas

On this screen, you can manage the company's areas. It offers the possibility of grouping users, making them belong to one sector.

■ Users

On this screen, you can manage the company's supervising users and operators. Bionova® Cloud will automatically inform by sending an email with a self-generated password to each new user created.

■ Connected incubators

All the incubators used by the system are listed here. In other words, all the incubators that were once connected to the machine.

It shows the connection status of the incubator. When you establish the connection to the incubator, it will be indicated in green.

■ Biological indicators

Consisting of two screens, it shows the data of the ongoing

and completed incubation processes. Furthermore, from this section, the user can complete all data to be manually uploaded; the table allows editing all relevant information about a BI result. Some of these fields will be automatically filled in with the information provided by the incubator, while others have to be filled in.

In process

This screen allows real-time monitoring of the incubation. It shows the time spent and the status of a biological indicator in the incubation process.

The user will find the list of fields (*columns*) below:

_Readout number: (*Automatic information*) Position of the incubator where the BI is being incubated and read.

_Position number: (*Automatic information*) Position of the incubator where the BI is being incubated and read.

_Program: (*Automatic information*) Selected program in the incubator.

_Start time: (*Automatic information*) shows the exact date and time of the incubation start.

_Time elapsed: (*Automatic information*) shows the minutes that have passed since the beginning of the incubation process. It should never be longer than the number of minutes indicated in the program.

_Incubator: (*Automatic information*) indicates the incubator that is performing the process.

Finished

The readouts that have completed the incubation process are displayed on the "Finished" screen. Now the readings can be filed. After the readouts are completed with the required information, it is possible to store the results.

You can view the filed results in the DASHBOARD > SAVED screen (see *Saved results*).

The colored icon at the beginning of each record represents the result of the readout.

Status colors:

Red: POSITIVE

Green: NEGATIVE

Grey: ANY OTHER STATUS

The user will find the list of fields (*columns*) below:

_Readout number: (*Automatic information*) Position of the incubator where the BI is being incubated and read.

_Position number: (*Automatic information*) Position of the incubator where the BI is being incubated and read.

↳ Biological indicators

- Start time: (*Automatic information*) Shows the exact date and time of the incubation start.
- End time: (*Automatic information*) Shows the exact date and time of the incubation end.
- Program: (*Automatic information*) Selected program in the incubator.
- Result: (*Automatic information*) Actual result given by the incubator (Positive, Negative, or Canceled).
- Average Temperature: (*Automatic information*) Average temperature during the BI incubation period.
- Ticket number: (*Automatic information*) Shows the device ticket number.
- Incubator: (*Automatic information*) indicates the incubator that is performing the process.

The screenshot shows a software interface for editing a reading. At the top, it says 'Edit Reading (109)'. Below that, there are several input fields:

- 'Product' dropdown: 'Choose a Product' with a 'Product List' button.
- 'Process' dropdown: 'Choose a Process' with a 'D Value' dropdown set to '1'.
- 'SCBI Condition' dropdown: 'Choose an SCBI Condition' with a 'Book Number' dropdown.
- 'Page Number' and 'Program Number' input fields.
- 'Cycle' dropdown: '1'.
- 'Sterilizer' dropdown: 'Choose a Sterilizer'.
- 'Term' and 'Notes' input fields.
- 'Load Number' input field.

At the bottom right, there is a green 'CONFIRM' button.

After the readouts are completed with the required information, it is possible to *File the results*.

Biological indicators for sterilization - Manual data

To edit a reading individually, see *Individual editing* point 7 | To edit multiple readings, see *Multiple editing*.

Product: BI being incubated.

Lot: BI lot number.

Process: conditions for the calculation of the D-value. The user can find this information in the product's "Instructions for Use" (which includes the BI lot quality certificate).

D-value: provides information about the resistance of that BI lot. The user can find this information in the product leaflet (which includes the BI lot quality certificate). Be careful to choose the correct D-value considering the "Process" selected in the previous field.

SCBI condition: You have to choose between 2 options: "Positive" (positive control, not exposed) or "Exposed" (to the sterilization cycle).

Book number: complete with the number of the record book corresponding to the sterilization process.

Page number: complete with the page number of the record book corresponding to the sterilization process.

Program number: the sterilizer program used for that cycle.

Cycle: cycle number, provided on the sterilizer display and/or the ticket.

Sterilizer: the sterilizer used for the process.

Conditions: you can fill in this box with details about the cycle (not mandatory information).

Comments: you can fill in this box with any relevant information not considered in the other fields (not mandatory information).

Load Number: (It can only be completed on the *Finished* screen) Day's load number.

Biological indicators for disinfection - Manual data

Product: BI being incubated.

Lot: BI lot number.

SCBI condition: You have to choose between 2 options: "Positive" (positive control, not exposed) or "Exposed" (to the sterilization cycle).

Book number: complete with the number of the record book corresponding to the sterilization process.

Page number: complete with the page number of the record book corresponding to the sterilization process.

Cycle: cycle number, provided on the sterilizer display and/or the ticket.

Disinfection machine: the disinfection machine used for the process.

Conditions: you can fill in this box with details about the cycle (not mandatory information).

Peroxide level:

Comments: you can fill in this box with any relevant information not considered in the other fields (not mandatory information).

The screenshot shows a software interface for editing a reading. At the top, it says 'Edit Reading (58)'. Below that, there are several input fields:

- 'Product' dropdown: 'Choose a Product' with a 'Product List' button.
- 'Process' dropdown: 'Choose a Process' with a 'Book Number' dropdown.
- 'SCBI Condition' dropdown: 'Choose an SCBI Condition' with a 'Desinfect' dropdown set to 'Choose a Desinfect'.
- 'Page Number' and 'Room ID' input fields.
- 'Term' and 'Notes' input fields.
- 'Room Volume' and 'Peroxide Concentration' input fields.

At the bottom right, there is a green 'CONFIRM' button.

Calculation of the sterilization process risk (SPR)

Bionova® Cloud allows calculating the sterilization process risk (SPR) of a set of results. You can access this function through the screen DASHBOARD > SPR in the side menu. This tool helps to make the most of the BI results, to control not only the sterility of a certain cycle but also the overtime performance of a specific sterilizer (called historical performance). This calculation is based on a very sophisticated and exclusive algorithm that takes into consideration the following information: cycle number, sterilizer (brand and serial number), number of BIs processed (in case you use more than one BI in a given cycle), the resistance of the BI lot (through the reported D-value), the BIs result (positive or negative) and the time to provide a positive result (in case you have a positive one), among others. The SPR calculation will be notified as None, Low, Medium, or High, following applicable international regulations. However, the user should consider that low risk does not mean it is not risky; a positive result implies that the controlled cycle cannot guarantee the load sterility, so the load needs to be

reprocessed, without exception. For this reason, the SPR value is not supposed to be examined cycle by cycle, but, on the contrary, considering a group of cycles over a considerable period. For this reason, these SPR calculations are used to represent the sterilizer performance using the SPR Graph tool. The user can perform a thorough inspection of each sterilizer at the CSSD by implementing this innovative tool, saving a lot of time for the sterilizer audit examination. In this way, the SPR provides useful data that can be used to reduce production costs, anticipate sterilizer repairs, or review sterilization procedures.

SPR Calculation Tab

When pressing the SPR Calculation tab, the program will display a form in which the user must indicate the grouping criteria for the calculation (sterilization cycle, sterilization machine). Once defined the parameters, the user must press the "Calculate SPR" button for the program to perform the calculation.

Note: Every institution must validate this tool according to the specific requirements for each implementation

The screenshot shows the RPE Chart interface with the 'Calculate RPE' tab selected. It includes fields for 'Sterilizer' (set to 'Sterilizer 1') and 'Cycle' (set to '1'). Below these, a green box displays the results: 'Reads analyzed: 3', 'Score: 67', and 'Score type: LOW'.


Historical Sterilizer Performance (SPR Graphic Tab)

By implementing this innovative tool, the user can perform a complete inspection and audit of each sterilizer at the CSSD. In this way, the SPR provides data that can be useful to

reduce production costs, anticipate sterilizer repairs, or review sterilization procedures.

The system places, by default, the current month in the date range. The sterilizer historical performance consists of a graph representing all SPR values for each cycle.


The screenshot shows the RPE Chart interface with the 'SPR Graphic' tab selected. It includes fields for 'Sterilizer' (set to 'Sterilizer 1'), 'Brand' (set to 'Ste 1'), and 'Serial' (set to '1111'). Below these, there are date range fields ('Date From: 01/01/2021' and 'Date To: 01/31/2021') and a 'GENERATE' button.

The X-axis represents the performance evolution through the cycle number, while the Y-axis draws the SPR score. Choose the sterilizer to be analyzed through the drop-down menu and press the "Generate" button. In this way, each sterilizer will have its historical performance.

Each cycle will have an SPR value represented by a point. The following are some examples of the evolution of the sterilizer performance and its interpretation to demonstrate the versatility and usefulness of this tool.


Example 1 | Correct the sterilizer performance.

In this case, the user can see that the sterilizer presented some fault, as shown in the first two SPR values with High levels (point 1). Thanks to this tool, it was possible to detect this sterilizer failure and to arrange proper maintenance.

The user can verify through the graph that after this maintenance, the sterilizer performance is acceptable since the SPR values return to Negative levels (point 2) (Cycle 3 to Cycle 9), confirming the correct operation from that moment.

Example 2 | Random sterilizer failure.


The sporadic occurrence of positive SPR values is an indicator of random failure, which could be due to premature cycle interruption, operator error, incorrect positioning of the BI within a test package or sterilization chamber, sterilizer chamber overload, among others. The detection of this type of failure could help improve the general protocol of the sterilization process through the analysis of the reasons for

these failures each time they occur.

If the user places the mouse pointer at any point, both the cycle number and the score received will be displayed, helping to identify the exact cycle. With this information and the results filtering tool, the user can track the rest of the data to evaluate the possible causes and reverse the situation.

Example 3 | Incorrect sterilizer functioning.

The sterilizer consecutively shows positive SPR levels, proving that there are no random situations but a consistent failure in some parameters.

In this case, the sterilizer showed normal SPR negative values and started giving positive values until it was repaired,

and the SPR values returned to negative levels. If you get a graph like this, you need to stop using this sterilizer urgently, schedule its maintenance/repair and re-evaluate its performance after repair (where you need to get consecutive negative SPR values).

■ Protein indicators

This screen shows all the saved results of the protein analysis, allowing you to sort and filter the results. It also shows the current PRO analysis, giving you the possibility to check the remaining time, to complete all the relevant information, and to save the result once finished.

In progress

This screen allows real-time monitoring of the incubation. It shows the elapsed time and the status of a protein indicator in the incubation process.

Below, the user will find the list of fields (*columns*).

_Reading number: (*Automatic information*) Position of the incubator where the BI is being incubated and read.

_Position number: (*Automatic information*) Position of the incubator where the BI is being incubated and read.

_Start time: (*Automatic information*) Indicates the exact date and time of the incubation start.

_Elapsed time: (*Automatic information*) Indicates the minutes that have passed since the incubation process began. It should never be longer than the number of minutes indicated in the program.

_Program: (*Automatic information*) Program selected in the incubator.

_Visual result: (*Manual information*) Visual confirmation of proteins through color change interpretation.

_Incubator: (*Automatic information*) Indicates which incubator is performing the process.

be filed. After the readouts are completed with the required information, it is possible to store the results.

You can view the filed results in the DASHBOARD > SAVED screen (see *Filed results*).

The colored icon at the beginning of each record represents the result of the readout.

Status colors:

Red: POSITIVE

Green: NEGATIVE

Grey: ANY OTHER STATUS

The user will find the list of fields (*columns*) below:

_Readout number: (*Automatic information*) Position of the incubator where the BI is being incubated and read.

_Position number: (*Automatic information*) Position of the incubator where the BI is being incubated and read.

_Start time: (*Automatic information*) Shows the exact date and time of the incubation start.

_End time: (*Automatic information*) Shows the exact date and time of the incubation end.

_Result: (*Automatic information*) Result of the protein test given by the incubator.

_Visual result: (*Manual information*) Visual confirmation of proteins through color change interpretation.

_Protein level (µg): (*Automatic information*) Quantitative Protein test result according to the incubator's information. This result is given in micrograms (µg) of protein. This information is provided through a BSA calibration curve.

_Ticket number: (*Automatic information*) Shows the device ticket number.

_Incubator: (*Automatic information*) Indicates the incubator that is performing the process.

Finished

The readouts that have completed the incubation process are displayed on the "Finished" screen. Now the readings can

↳ Protein indicators

Information to be completed manually

To edit a reading individually, see *Individual editing* point 7|

To edit multiple readings, see *Multiple editing*.

The following is the form for completing manual information:

Edit Protein Indicator Reading (90)

Product Choose a Product	Product lot _____
Book Number _____	Page Number _____
Cycle _____	Surface _____
Washer Choose a Washer	Visual Result _____
Observations _____	
CONFIRM	

_Product: BI being incubated.

_Lot: BI lot number.

_Book number: complete with the number of the record book corresponding to the sterilization process.

_Page number: complete with the page number of the record book corresponding to the sterilization process.

_Cycle: cycle number, provided on the sterilizer display and/or the ticket.

_Surface: provides information about the item that has been sampled.

_Washing machine: the machine used for the process.

_Visual result: Visual confirmation of proteins through color change interpretation.

_Comments: you can fill in this box with any relevant information not considered in the other fields (not mandatory information).

Development of the baseline

The first step concerns the sampling of reprocessed surgical instruments representing the full range of your workload to provide the basis on which the monitoring system can be developed (benchmark).

There are two main approaches to instrument sampling:

_If a SINGLE measurement takes place at any time

_Measurements that are naturally divided into GROUPS, established by the user (for example, five instruments per week).

In the case of the SINGLE approach, a total of 20 measurements should initially be sufficient to estimate the baseline. In the case of GROUPS, it will take more than 30 measurements.

Procedure

1.1| Process the instrument in your washing machine as usual, following the recommendations of both the instrument and the washing machine manufacturers.

1.2| Select a representative instrument to test for residual protein.

1.3| Test the selected instrument with the Chemodye® PRO1MICRO device, strictly following the device's instructions for use.

1.4| Incubate the PRO1MICRO used, according to the indicator.

1.5| Before or during incubation, connect the incubator to your PC following the instructions detailed in the incubator user manual. Make sure you have the Bionova Cloud® Agent installed on your PC.

1.6| Once inside the Traceability Software website, go to the HYGIENE PROTEIN SYSTEM > IN PROGRESS screen.

1.7| The type of monitoring and results are automatically completed, as well as the protein level and limit. The visual result, however, should be completed according to the color of the reactive solution after the incubation.

1.8| The next step is to complete all the manual information needed (see *Information to be completed manually*).

Caution: To save the visual result, remove the protein pen, and check the solution color as soon as the quantitative reading is completed.

The readings can now be "filed." After completing the readings with the required information, it is possible to *File the results*.

HTM 01-01

Controlling the reprocessing of instruments through statistical control charts assures that the process is working as anticipated, and the residual proteins in the tools are within the expected range.

You can implement an Internal Quality Assurance System (IQAS) to control, over time, the effectiveness of the washing process.

Control analysis

Continuous monitoring (Recommendation HTM 01-01):

You should create a graph for each Chamber/Washer-Disinfector. The number and type of instruments should be tested as follows:

50 instruments per Chamber/Washer-Disinfector, at least every three months, chosen among those that are difficult to clean (for example, box joints, serrated edges, hinges, graters, and complex surfaces) when used. Other difficult to clean instruments should also be identified and included in this test.

Procedure

1.1 | Once the first 20 measurements have been obtained, the control analysis can be started by going to the screen DASHBOARD > HTM 01-01

The following form will be displayed on the screen:

The screenshot shows a dropdown menu for 'Machine' set to 'WASHA'. Below it are two dropdown menus for 'Date range': 'From 01/01/2021' and 'Until 01/31/2021'. Under 'Measurement Grouping', there is a dropdown set to '1' and another set to 'Individual Chart'. A green 'GENERATE' button is at the bottom.

1.2 | Select your washing machine by clicking on the drop-down menu. By default, the date range for the analysis corresponds to the current month.

1.3 | Modify the measurement grouping option if you wish to execute a group approach. According to HTM 01-01, more than 30 measurements should be considered for the basis, and 10 as the maximum number per group.


1.4 | If you wish to continue with individual focusing (measurement grouping = 1), two options will be displayed in the **Graphics** menu:

— "Graph I" (Individual Graph): log10 protein residues vs. measurement number.

— "Mr Graph" (Moving Range Graph): absolute differences between log10 protein residues vs. the number of differences.

1.5 | If you want to execute a group approach (measurement grouping >01) two options will be displayed in the **Graphics** menu:

Example 1 | Out of control type 1: Instability. In this case, 8 or more consecutive measurements are on the same side of the average.

— "XBAR Graph": group average vs. group number.

— "R-Graph" (Range graph): group range vs. group range number.

Results interpretation

The monitoring of each process should vary randomly within the average and control limits (Upper Control Limit and Lower Control Limit) and the process variation should also be random and within the control limit. If a measurement exceeds the Upper Control Limit, the sample must be washed again and the instruments sampled again. The maximum upper limit allowed by the standard is 5 µg for each side of the instrument (according to HTM 01-01).

However, measurements that are outside the control limits are an indication that the process is no longer under control. This may warrant an investigation to detect and eliminate any underlying cause.


Before interpreting Graph I (individual approach) or XBAR Graph (group approach), the MR Graph or R Graph (respectively) should be examined to determine if the process variations are under control. If the MR-chart or the R-chart are not controlled, then the control limits of Graph I or XBAR cannot be considered accurate. Even when the measurements of the Graph I or XBAR are within the control limits, it does not necessarily mean that the system is under control.

It is important to analyze any systematic behavior in sequential measurements. For example, if there were 10 successful measurements, all on the centerline, then this may indicate that has occurred a systematic change in the process, as we would expect half of these measurements to be divided on each side of the centerline. To assist in the interpretation of statistical control process graphs, several situations have been suggested (HTM 01-01) that may indicate an "out of control" process. Our software is prepared to monitor the "out of control" programs, detecting these situations indicated below:

Example 2| Out of control type 2: Alternating. In this case, there are 14 measurements where there is an alternating pattern.

Example 3| Out of control type 3: Monotonous. In this case, 6 consecutive measurements show a monotonous trend (all rising or falling).

■ Chemical indicators

This tab allows the detection and filing of chemical indicators of washing and sterilization through the use of the Automatic System for Quality Control and Traceability of Washing,

Disinfection, and Sterilization Processes Trazanto®. In this section, you can also consult previously filed results. It is composed of two screens that differentiate the readings of sterilization indicators and washing machines.

Detection

CHEMICAL INDICATORS COMPATIBLE WITH TRAZANTO		1																																																					
STERILIZER PROVISOARY READINGS																																																							
<input type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>		2 SCAN																																																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>READ NUMBER</th> <th>TEST DATE TIME</th> <th>RESULT</th> <th>MANUFACTURE DATE</th> <th>VISUAL RESULT</th> <th>PRODUCT</th> <th>PRODUCT LOT</th> <th>BRAND</th> </tr> </thead> <tbody> <tr> <td>7</td> <td>10/16/2020 14:19:42</td> <td>REJECTED</td> <td>07/03/2019</td> <td></td> <td>CDWA4JAY</td> <td>CDWA4030719</td> <td>CHEMDYE</td> </tr> <tr> <td>6</td> <td>10/16/2020 14:13:48</td> <td>REJECTED</td> <td>07/03/2019</td> <td></td> <td>CDWA4JAY</td> <td>CDWA4030719</td> <td>CHEMDYE</td> </tr> <tr> <td>5</td> <td>10/14/2020 07:34:45</td> <td>ACCEPTED</td> <td>02/05/2020</td> <td></td> <td>CDWE</td> <td>IT261YS050220</td> <td></td> </tr> <tr> <td>4</td> <td>10/14/2020 07:30:24</td> <td>ACCEPTED</td> <td>02/05/2020</td> <td></td> <td>CDWE</td> <td>IT261YS050220</td> <td></td> </tr> <tr> <td>3</td> <td>10/14/2020 07:27:56</td> <td>ACCEPTED</td> <td>02/05/2020</td> <td></td> <td>CDWE</td> <td>IT261YS050220</td> <td></td> </tr> </tbody> </table>								READ NUMBER	TEST DATE TIME	RESULT	MANUFACTURE DATE	VISUAL RESULT	PRODUCT	PRODUCT LOT	BRAND	7	10/16/2020 14:19:42	REJECTED	07/03/2019		CDWA4JAY	CDWA4030719	CHEMDYE	6	10/16/2020 14:13:48	REJECTED	07/03/2019		CDWA4JAY	CDWA4030719	CHEMDYE	5	10/14/2020 07:34:45	ACCEPTED	02/05/2020		CDWE	IT261YS050220		4	10/14/2020 07:30:24	ACCEPTED	02/05/2020		CDWE	IT261YS050220		3	10/14/2020 07:27:56	ACCEPTED	02/05/2020		CDWE	IT261YS050220	
READ NUMBER	TEST DATE TIME	RESULT	MANUFACTURE DATE	VISUAL RESULT	PRODUCT	PRODUCT LOT	BRAND																																																
7	10/16/2020 14:19:42	REJECTED	07/03/2019		CDWA4JAY	CDWA4030719	CHEMDYE																																																
6	10/16/2020 14:13:48	REJECTED	07/03/2019		CDWA4JAY	CDWA4030719	CHEMDYE																																																
5	10/14/2020 07:34:45	ACCEPTED	02/05/2020		CDWE	IT261YS050220																																																	
4	10/14/2020 07:30:24	ACCEPTED	02/05/2020		CDWE	IT261YS050220																																																	
3	10/14/2020 07:27:56	ACCEPTED	02/05/2020		CDWE	IT261YS050220																																																	
ROWS PER PAGE 5 < 1-5/6 > >																																																							

↳ Chemical indicators

- 1 | The drop-down menu, when clicked, specifies which indicators are compatible with Trazanto®.
- 2 | Press the  Analyze button to automatically detect and read an indicator. A detection screen like the following will be displayed:

Important: Make sure you have the Bionova® Cloud Agent installed so that your Trazanto reader can connect to the Bionova® Cloud web.

Before using the tool, the Trazanto® reader must be turned on and connected to the computer through the USB cable provided with it, and the chemical indicator to be checked should be placed inside the reader (maximum 1 at a time). For further details on how to install and use this device, please refer to the manual that came with it. This screen will inform you of any problems detected.

Once the detection and reading is finished, you can enter the Bionova® Cloud website, displaying the reading you just made on the corresponding **Chemical indicators** screen.

Information to be completed manually

To edit a reading individually, see *Individual editing* point 7 | To edit multiple readings, see *Multiple editing*.

Cycle number: cycle number of the sterilization or washing process.

Temperature: temperature of the sterilization or washing process cycle.

Machine [Sterilizer or Washer]: descriptive name of the

sterilizer or washer used in the process.

Machine serial number: serial number of the sterilizer or washing machine. This field is completed by the machine search tool .

Exposure time: time during which the indicator was exposed to the process.

Package number (Exclusive for Sterilization processes): The number of the package in which the chemical indicator is found at the time of sterilization.

Location (Exclusive for washing processes): Location of the chemical indicator inside the washing machine.

Relative dampness (Exclusive for Sterilization processes): Relative humidity condition, given in percentage, in which the sterilization process was performed.

Concentration: percentage of concentration of the detergent (for washing) or the sterilizing agent (for sterilization processes).

Type of detergent (Exclusive for washing processes): You can select between enzymatic, alkaline, non-ionic, or you can complete with the value you want.

Program: program number of the equipment used for that cycle.

Visual result: visual confirmation of the process result.

Book Number: number of the record book corresponding to the process.

Page number: number of record page corresponding to the process.

Observations: you can fill in this box with any relevant information that has not been considered in the other fields. This field can be filled in or modified after the result has been saved.

Water hardness (Exclusive for Washing processes): Parameter that defines the quality of the water. This data can be obtained from the quality parameters of the water treatment system used in the establishment.

Water pressure (Exclusive for Washing processes): Pressure of the water used in the cycle. You can obtain this data from the washing equipment.

After the readouts are completed with the required information, it is possible to *File the results*.

You can view the filed results in the DASHBOARD > SAVED screen (see *Filed results*).

■ Equipment

In these screens, you can control the sterilizers, washing machines, and disinfectors' registrations, deregistrations, and

modifications. They are the pieces of equipment that belong to a company, and later it will be indicated in a reading which equipment carried out the process.

■ Licenses

My licenses

An active license is required to use certain Bionova® Cloud modules or screens. The Administrator user (Admin) has the permission to request a license.

To know in more detail the structure and internal format of these files for a correct implementation, it is available the document: Manual for the developer.

Menu item: LICENSE > MY LICENSES

You will be able to view the status of your licenses, start and end date.

There are three licenses to purchase:

1|Exporter: used to export reading records in "Filed" status.

2|Bionova® Exchange: needed to subscribe to the Bionova® Cloud readings notification service.

3|Bionova® Exchange+(Plus): same as a normal subscription. If it is the first time you will request a license, the system will activate the Trial automatically.

■ Filed results

The DASHBOARD > SAVED screen contains a table showing all the information that was completed during the indicator's incubation and saved when this incubation concluded. Only the BI readings can be edited, changing the visual result.

In the upper-right side of the screen, you can select the type of indicator you want to consult.

i|Biological Indicators

- (a) Sterilization
- (b) Disinfection

ii|Protein indicators

iii|Chemical Indicators

- (a) Cleaning
- (b) Sterilization

Biological Indicator ▾ Protein Indicator Chemical Indicators ▾

The system performs a reading search, differentiating the type of user: administrator and supervisor.

When the user is an Administrator:

You can view all of the company's indicator readings.

When the user is a Supervisor:

Only will be visible those readings of indicators from the sector to which the user belongs.

Filter to search for specific information on historical readings

■ Reports download

By default, the download includes all Saved results. In this case, Bionova® Cloud will notify you that the export may take a long time. However, if a filter was previously applied, the export will be restricted to the records that match the filtering criteria. The system will download a pdf file containing the results when you click on the following button: . You could later print the report from the browser.

■ Equipment calibration

Before starting the calibration process, make sure that you will not need the incubator for the next 2-3 hours because that is approximately how long this process takes.

The temperature calibration tool available in the Bionova® Cloud Agent program allows you to adjust the temperature of the Bionova® incubator. It is an easy-to-use tool that facilitates regular temperature calibration according to a procedure determined by the program. As a consequence of the aging of the equipment, some changes occur in the heating system. To compensate for these small changes, this function allows a precise adjustment of the incubation temperature. This new tool will allow generating a calibration certificate for the unit. Some information recorded in this certificate makes reference to the institution to which the equipment belongs, the staff who

← **Equipment calibration**

performed the procedure, the information about the equipment, the identification of the reference device, the environmental conditions, and the calibration data.

A failure in the calibration process or operating errors could cause the equipment to malfunction. Be sure you have the knowledge and equipment to perform this procedure safely. Terragene S.A. will not be responsible for any loss or damage that may result from improper use of the equipment or failure of the calibration tool.

Calibration process

To perform this procedure, make sure you have an external thermometer suitable for the thermometer hole located in the back of the incubator. Our digital thermometer TB-IC1020 has been developed with the appropriate specifications.

In case you use a different brand, make sure that the specifications match those of TB-IC1020.

To start the calibration process, press the button  and follow the instructions in the program. Start the calibration process only when the temperature of the equipment is stable.

Once the calibration is completed, a calibration certificate will automatically appear on the screen, and you can print or save it.

■ **Other functions**

Clock synchronization

By pressing the icon  you can synchronize the incubator clock (*date and time*) with your computer system clock.

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