Bionova® Hyper Auto-reader

Symbols | Símbolos | Símbolos

You may see one or more of these symbols on the packaging or labeling of this product. | Es posible que vea uno o más de estos símbolos en el empaque o el etiquetado de este producto. | Poderá ver um ou mais destes símbolos na embalagem ou rotulagem deste produto.

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EN Bionova® Hyper Auto-reader

Product description

The Bionova® Hyper Auto-reader was designed for the incubation and automatic readout of the Bionova® Hyper Self-contained Biological Indicator (SCBI) for monitoring plasma or vaporized hydrogen Peroxide (H2O2) sterilization processes. The Bionova® Hyper Auto-reader has two independent readout positions. Each position automatically detects when an SCBI has been placed for incubation and delivers a readout result after a few minutes using advanced fluorescence techniques. The Bionova® Hyper Auto-reader detects positive and negative SCBIs easily and quickly. A positive result can also be shown by a change in the culture medium color when extended incubations are carried out. Read the Instructions for Use of the SCBI for more information. The option to perform or not an extended incubation depends on the internal protocols of each laboratory or hospital. The Bionova® Hyper Auto-reader has built-in USB, Wi-Fi, and Bluetooth® technology connectivity. The electronic tickets (e-tickets) created for each readout process can be accessed using a compatible device for easy documentation compliance and traceability.

Indications for Use

Terragene S.A. recommends the Hyper Auto-reader Incubator (BHY) incubates at 60 °C and reads the Terragene® Bionova® Hyper SCBIs at the times prescribed in the User Manual.

Safety Information

-  To avoid risks and/or damaging the device:
 - Use indoors only.
 - Do not place the auto-reader in a room exposed to direct sunlight or to high luminous intensity lamps.
 - Do not place the auto-reader near devices that emit strong electromagnetic fields.
 - Do not use the auto-reader on sloping surfaces or on surfaces that have been knocked, that may vibrate, be exposed high temperatures or high relative humidity.
 - Disconnect the power cord before cleaning.
 - Do not use abrasive, corrosive cleaners or disinfectants.
 - Do not immerse into any liquid. Do not pour liquid into it.
 - Make sure the auto-reader is connected to an appropriate electrical mains outlet socket.
 - Use only the included power supply (AC power adapter), power supply's AC plug, power supply cords, and USB cable. Check that all the items included are in good condition every day. If any of them is damaged, stop using it. The use of cables, adapters, cords, and/or power supplies different from the ones included may cause fires, electrical shocks, or even physical injuries.
 - Do not plug any devices into the auto-reader's USB port other than a personal computer (PC). The PC must be compliant with IEC 60950-1, IEC 62368-1 or comparable, with safety extra-low voltages on its USB ports. Ask a qualified technician to check that the device is compatible. Attaching any other device to the USB port may damage the auto-reader and may not be safe for the user.
 - Do not attempt to repair the auto-reader yourself. This could lead to major and irreversible damages to the device. If the auto-reader malfunctions, contact your local distributor for further assistance.

-  To reduce the risk of using incompletely sterilized loads:
 - Please read, understand, and follow the Instructions for Use of each SCBI before incubating it.
 - Do not remove the SCBI before the incubator shows the final readout result.
 - Check that the spore carrier has been completely wetted by the culture medium.

-  To avoid the risk of injury by glass fragments produced when crushing the glass ampoule inside the SCBI tube:
 - Cool the SCBI during the indicated time before crushing the ampoule.
 - Do not handle the SCBI excessively since this might cause the glass ampoule to explode.
 - Wear safety gloves and goggles while removing the SCBI from the sterilizer, pressing the SCBI's cap, and crushing the SCBI's ampoule.
 - Do not use your finger to crush the SCBI. Use the ampoule crusher instead.

-  To avoid a potentially hazardous situation:
 - Avoid contact block inside each incubation position.
 - Do not insert your fingers, or any other element, into the incubation positions.
 - Place only compatible indicators inside the incubation positions.

-  To avoid SCBIs from absorbing fluorescent particles:
 - Avoid direct contact between the SCBIs and chemical indicators or tapes before the SCBI incubation.
 - Terragene S.A. recommends the Hyper Auto-reader should not be handled by the indicator.

Note: Do not use this product other than in the manner specified by Terragene S.A., otherwise the protection provided by the product might be affected.
Note: Only Terragene S.A. authorized personnel may access or service the internal components of the auto-reader. Parts or components inside the auto-reader should not be handled by the user.

Operating conditions

Input parameters	Operating conditions	Units
Voltage range	(100-240)	AC Volts
Frequency	50-60	Hertz
Current	0.6	Amperes
Output parameters	Values	Units
Voltage	12	DC Volts
Current	1.5	Amperes

Terragene S.A. recommends the use of Uninterruptable Power Supply (UPS) instead of voltage stabilizers, since they fulfill two functions: to stabilize and maintain the energy during a power outage.

Environmental operating conditions	Operating conditions	Units
Altitude	3500 (max.)	Meters
Operation temperature	10-30	Celsius
Relative humidity	30-80	%
Installation/Over-voltage	Category II	-
Pollution degree	2	-
Storage temperature	5-40	Celsius
Voltage	12	DC Volts

Device safety compliance

The auto-reader complies with the following standards and directives:

Electrical Safety	IEC 61010-1 IEC 61010-2-010 Low Voltage Directive 2014/35/EU
Electromagnetic Immunity	IEC 60601-1-2 EN 60601-1-2
Electromagnetic Compatibility (EMC)	EN 61326-1 EN 62311 RED Directive 2014/53/EU 47 CFR Part 15
European Commission	RoHS Directive 2011/65/EU WEEE Directive 2012/19/EU

EU Wireless Compliance: The auto-reader complies with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following Internet address: www.terragene.com/DOC

RF transmitter specifications: The Wi-Fi transmitter total output power is 19.9 dBm E.I.R.P (97.72 mW) in the frequency band of 2.4 GHz ISM band / 2.412 to 2.462 GHz (Channels 1 to 11). The Bluetooth® transmitter total output power is 4.6 dBm E.I.R.P (2.88 mW) in the frequency band of 2.4 GHz ISM band / 2.402 to 2.480 GHz.

FCC/IC Wireless Compliance: The auto-reader complies with applicable FCC/IC directives.

It contains a certified transmitter module: FCC ID: 2AC7Z-ESPWROOM32D / IC ID: 21098- ESPWROOM32D.

This auto-reader complies with Part 15, Subpart A of the FCC rules: 15.247 Operation within the bands 902- 928 MHz, 2400-2483.5 MHz and 5725-5850 MHz, and 15.209 Radiated emission limits.

General requirements, and the Canadian standards: RSS-Gen General Requirements for Compliance of Radio Apparatus, and RSS-247 Digital Transmission Systems (DTSS), Frequency Hopping Systems (FHSs) and License-Exempt Local Area Network (LE-LAN) Devices. It may only be used if the following two conditions are met:
1. This auto-reader may not cause interference.
2. This auto-reader must tolerate interference, including interference that may cause undesired operation of the device.

The auto-reader does not represent a photobiological risk and does not generate dangerous optical radiation if used within normal operation conditions as per the requirements of IEC 62471 Standard.

Designed under Quality Management System standards ISO 13485:2016/EN ISO 13485:2016.

For more information, please refer to the User Manual.

Note: The auto-reader has been evaluated for conformity for use in business environments. There may be concerns about radio interference if used in home environments. This applies exclusively to the Korean market due to KC certification requirements.

Cleaning and maintenance

Cleaning and decontamination of external surfaces

Disconnect the power supply cable and the USB cable from the auto-reader. If the device is hot, wait until it has cooled down before handling it.
Clean the external surfaces of the device using a microfiber cloth moistened with a solution of mild dish washing detergent and water. Wring out the cloth so it is damp but not dripping before cleaning and wipe the outer surfaces of the auto-reader. Then, moisten a clean microfiber cloth with just water, and repeat the procedure until all traces of detergent are removed from the external surfaces of the device. After cleaning, allow the auto-reader to air dry for at least 1 hour before connecting the power supply cable or the USB cable again.
This cleaning procedure can be followed whenever considered appropriate.
This cleaning procedure must be followed each time a spill occurs on any of the external surfaces of the auto-reader.
If further cleaning is required, or if you have doubts about the cleaning agents you may use, please contact your local distributor.

Note: Do not clean the internal parts of the device.

Note: Do not pour or liquid on o the device or immerse it in any liquid. Do not allow any liquid to enter the device while it is being cleaned.

Maintenance

The auto-reader does not require routine maintenance.

Warranty

The products are guaranteed to be free from material and workmanship defects when properly installed, maintained, and used for the purpose indicated in the applicable product label and/or the User Manual. The warranty only applies to the original purchaser

Term

The warranty period for the Bionova® Hyper Auto-reader is 1 (one) year from the date installation and may never be extended beyond 5 (five) years from the product's date of manufacture.

Limitation of liability

Terragene S.A. shall not be held liable for any loss or damage that result from the unsuitable use of the equipment, negligence or user's full responsibility.

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EN Instructions for use

Note: Remember to set-up the auto-reader local date and time before starting incubations.

Setting Up your auto-reader

Place the auto-reader on a firm surface, free from vibrations, away from direct sunlight, hot or cold air flows, chemical and corrosive or flammable substances. Do not place the auto-reader in a way that would make it difficult to remove supply's AC plug from mains. Leave a space of at least 10 cm from the auto-reader to the closest wall. Do not move the auto-reader periodically while it is being used. Connect the auto-reader to a secure and stable electrical mains outlet socket. Do not wet or heat the auto-reader. If liquid is spilled on the device, disconnect it and dry it immediately. Before turning it on, check that all reading positions are empty.

Using your auto-reader

1| Switch on your auto-reader by connecting the power supply's AC plug to the mains and then connect the plug at the other end of power supply to the rear of the auto-reader. (See image 1)
2| Set the auto-reader local date and time running the Bionova® Wireless Assistant Application on a compatible mobile device.
3| Wait until the Temperature & stability indicator light stops blinking. Place an SCBI in a readout position. The incubation will start automatically. (See image 2)
4| Wait until the incubation has finished. The position status indicator light will turn green to indicate a negative result, or red to indicate a positive result. If a positive result is obtained, press the Configuration & alarm cancellation button to silence the alarm. Do not remove the SCBI until a result has been produced. (See image 3)
5| Run Bionova® Cloud Web Application or the Bionova® Wireless Assistant Application to check the incubation result e-ticket.
Disposal: Dispose of the SCBIs in accordance with your country's waste-disposal guidelines.

Connectivity

To connect your auto-reader to a personal computer (PC) by USB, connect one end of the USB cable provided to the rear of the auto-reader and the other end to an USB port on your PC. Connecting the auto-reader to a PC to use Bionova® Cloud Web Application. For more information, please refer to the User Manual.
To access the extended features of your auto-reader using Wi-Fi or Bluetooth®, please download the Bionova® Wireless Assistant Application to your preferred mobile device and follow the steps detailed on the Application's tutorial.
Note: This auto-reader is only compatible with 2.4 GHz Wi-Fi networks.
Note: For cybersecurity reasons, the auto-reader can only be connected to secure, password-protected Wi-Fi networks. It cannot be connected to open or public networks.

For a detailed overview on the connection process, please refer to the User Manual.

Connectivity credentials

During the connection process, the following credentials may be useful:

Device Internal Wi-Fi network credentials
Wi-Fi network name: HYPER_[batch][serial number]
Wi-Fi network password: abcd1234

Bluetooth® credentials
Device Bluetooth® ID: HYPER_[batch][serial number]

Product compatibility and incubation programs

Auto-reader incubation programs

The following incubation program is available for your auto-reader.

Temperature selection	Time selection	Incubation program
60 °C	Hyper (5 min.)	Hyper at 60 °C

Note: Please read and follow the SCBI Instructions for Use before using it.

Troubleshooting

Refer to the User Manual for detailed troubleshooting measures
Please visit www.terrakene.com/BIONOVA-HYPER or scan the QR code to access the auto-reader User Manual.

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