Bionova® Photon Auto-reader

Symbols | Símbolos | Símbolos

You may see one or more of these symbols on the packaging or labeling of this product. | Es posible que vea uno o más de estos símbolos en el emvase o el etiquetado de este producto. | Poderá ver um ou mais destes símbolos na embalagem ou rotulagem deste produto.

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EN Bionova® Photon Auto-reader

Product description

The Bionova® Photon Auto-reader was designed for the incubation and automatic readout of the Bionova® Photon Self-contained Biological Indicator (SCBI) for monitoring steam sterilization cycles. The Bionova® Photon Auto-reader has two independent readout positions. Each position automatically detects when an SCBI has been placed for incubation and delivers a readout result after a few seconds using advanced fluorescence techniques. The Bionova® Photon Auto-reader detects positive and negative SCBIs easily and quickly. A positive result can also be shown by a change in the culture medium color when extended incubations are carried out. Read the Instructions for Use of the SCBI for more information. The option to perform or not an extended incubation depends on the internal protocols of each laboratory or hospital. The Bionova® Photon Auto-reader has built-in USB, Wi-Fi, and Bluetooth® technology connectivity. The electronic tickets (e-tickets) created for each readout process can be accessed using a compatible device for easy documentation compliance and traceability. The Bionova® Photon Auto-reader is intended to be marketed as an Over-the-Counter (OTC) device.

Indications for Use

The Bionova® Photon Auto-reader Incubator (BPH) incubates at 60 °C and reads the Terragene Bionova® Photon SCBIs at the times prescribed in the User Manual.

Safety Information

To avoid risks and/or damaging the device:

- Use indoors only.
- Do not place the auto-reader in a room exposed to direct sunlight or to high luminous intensity lamps.
- Do not place the auto-reader near devices that emit strong electromagnetic fields.
- Do not use the auto-reader on sloping surfaces or on surfaces that might be knocked, that may vibrate, be exposed to high temperatures or high relative humidity.
- Disconnect the power cord before cleaning.
- Do not use abrasive, corrosive cleaners or disinfectants.
- Do not immerse into the auto-reader. Then, moisten a clean microfiber cloth with water only, and repeat the procedure until all traces of detergent are removed from the external surfaces of the device. After cleaning, allow the auto-reader to air dry for at least 1 hour before connecting the power supply cable or the USB cable again.
- This cleaning procedure can be followed whenever considered appropriate.
- This cleaning procedure must be followed each time a spill occurs on any of the external surfaces of the auto-reader.
- If further cleaning is required, or if you have doubts about the cleaning agents you may use, please contact your local distributor.
- Note:** Do not clean the internal parts of the device.
- Note:** Do not pour liquid on the device or immerse it in any liquid. Do not allow any liquid to enter the device while it is being cleaned.

To reduce the risk of using incompletely sterilized loads:

- Please read, understand, and follow the Instructions for Use of each SCBI before incubating it.
- Do not remove the SCBI before the auto-reader shows the final readout result.
- Check that the spore carrier has been completely wetted by the culture medium.

To avoid the risk of injury by glass fragments produced when crushing the glass ampoule inside the SCBI tube:

- Cool the SCBI during the indicated time before crushing the ampoule.
- Do not handle the SCBI excessively since this might cause the glass ampoule to explode.
- Wear safety gloves and goggles when removing the SCBI from the sterilizer, pressing the SCBI's cap, and crushing the SCBI's ampoule.
- Do not use your finger to crush the SCBI. Use the ampoule crusher instead.

To avoid a potentially hazardous situation:

- Avoid contact with the hot metal block inside each incubation position.
- Do not insert your fingers, or any other element, inside the incubation positions.
- Place only compatible indicators inside the incubation positions.

To avoid SCBIs from absorbing fluorescent particles:

- Avoid direct contact between the SCBIs and chemical indicators or tapes before the SCBI incubation.
- Avoid excessive SCBI handling that may lead to fingerprints or glove talc imprinted over the indicator.

Note: Do not use this product other than in the manner specified by Terragene S.A., otherwise the protection provided by the product might be affected.

Note: Only Terragene S.A. authorized personnel may access or service the internal components of the auto-reader. Parts or components inside the auto-reader should not be handled by the user.

Operating conditions

Input parameters	Operating conditions	Units
Voltage range	(100-240)	AC Volts
Frequency	50-60	Hertz
Current	0.6	Amperes
Output parameters	Values	Units
Voltage	12	DC Volts
Current	1.5	Amperes

Terragene S.A. recommends the use of Uninterruptable Power Supply (UPS) instead of voltage stabilizers since they fulfill two functions: to stabilize and maintain the energy during a power outage.

Environmental operating conditions

Environmental conditions	Operating conditions	Units
Altitude	3500 (max.)	Meters
Operation temperature	10-30	Celsius
Relative humidity	30-80	%
Installation/Over-voltage	Category II	-
Pollution degree	2	-
Storage temperature	5-40	Celsius
Voltage	12	DC Volts

Device safety compliance

The auto-reader complies with the following standards and directives:

Electrical Safety	IEC 61010-1 IEC 61010-2-010 Low Voltage Directive 2014/35/EU
Electromagnetic Immunity (EMI)	IEC 60601-1-2 EN 60601-1-2 EN 61326-1 RED Directive 2014/53/EU 47 CFR Part 15
Electromagnetic Compatibility (EMC)	EN 61326-1 EN 62311 RED Directive 2014/53/EU 47 CFR Part 15
European Commission	RoHS Directive 2011/65/EU WEEE Directive 2012/19/EU

EU Wireless Compliance: The auto-reader complies with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following Internet address: www.terragene.com/DOC

RF transmitter specifications: The Wi-Fi transmitter total output power is 19.9 dBm E.I.R.P (97.72 mW) in the frequency band of 2.4 GHz ISM band / 2.412 to 2.462 GHz (Channels 1 to 11). The Bluetooth® transmitter total output power is 4.6 dBm E.I.R.P (2.88 mW) in the frequency band of 2.4 GHz ISM band / 2.402 to 2.480 GHz.

FCC/IC Wireless Compliance: The auto-reader complies with applicable FCC/IC directives. It contains a certified transmitter module: FCC ID: 2AC7Z-ESPWROOM32D / IC ID: 21098- ESPWROOM32D.

This device complies with Part 15, Subpart A of the FCC rules: 15.247 Operation within the bands 902- 928 MHz, 2400-2483.5 MHz and 5725-5850 MHz, and 15.209 Radiated emission limits.

General requirements and the Canadian standards: RSS-Gen General Requirements for Compliance of Radio Apparatus, and RSS-247 Digital Transmission Systems (DTSS), Frequency Hopping Systems (FHSS) and License-Exempt Local Area Network (LE-LAN) Devices. It may only be used if the following two conditions are met:
1. This auto-reader may not cause interference.
2. This auto-reader must tolerate interference, including interference that may cause undesired operation of the device.

The auto-reader does not represent photobiological risk and does not generate dangerous optical radiation if used within normal operation conditions as per the requirements of IEC 62471 Standard.

Designed under Quality Management System standards ISO 13485:2016 / EN ISO 13485:2016. For more information, please refer to the User Manual.
Note: The auto-reader has been evaluated for conformity for use in business environments. There may be concerns about radio interference if used in home environments. This applies exclusively to the Korean market due to KC certification requirements.

Cleaning and maintenance

Cleaning and decontamination of external surfaces

Disconnect the power supply cable and the USB cable from the auto-reader. If the device is hot, wait until it has cooled down before handling it. Clean the external surfaces of the device using a microfibr cloth moistened with a solution of mild dish washing detergent and water. Wring out the cloth so it is damp but not dripping before cleaning and wipe the outer surfaces of the auto-reader. Then, moisten a clean microfiber cloth with water only, and repeat the procedure until all traces of detergent are removed from the external surfaces of the device. After cleaning, allow the auto-reader to air dry for at least 1 hour before connecting the power supply cable or the USB cable again. This cleaning procedure can be followed whenever considered appropriate. This cleaning procedure must be followed each time a spill occurs on any of the external surfaces of the auto-reader.

If further cleaning is required, or if you have doubts about the cleaning agents you may use, please contact your local distributor.

Note: Do not clean the internal parts of the device.
Note: Do not pour liquid on the device or immerse it in any liquid. Do not allow any liquid to enter the device while it is being cleaned.

Maintenance

The auto-reader does not require routine maintenance.

Warranty

The products are guaranteed to be free from material and workmanship defects when properly installed, maintained and used for their intended purpose indicated in the applicable product label and/or the User Manual. The warranty only applies to the original purchaser.

Term

The warranty period for the Bionova® Photon Auto-reader is 1 (one) year from the date of installation and may never be extended beyond 5 (five) years from the products' date of manufacture.

Limitation of liability

Terragene S.A. shall not be held liable for any loss or damage that result from the unsuitable use of the equipment, negligence or user's full responsibility.

Technical assistance

For information regarding the functioning or condition of the product, contact Terragene S.A. directly to info@terragene.com or through our website www.terragene.com. In USA, please call +1 844 837 7243. Terragene S.A. Ruta Nacional N° 9, Km 280 - CP 2130, Parque Industrial Micropi-Alvear-Santa Fe-Argentina.

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EN Instructions for use

Note: Remember to set-up the auto-reader local date and time before starting incubations.

Setting Up your auto-reader

Place the auto-reader on a firm surface, free from vibrations, away from direct sunlight, hot or cold air flows, chemical and corrosive or flammable substances. Do not place the auto-reader in a way that would make it difficult to remove supply's AC plug from mains. Leave a space of at least 10 cm from the auto-reader to the closest wall. Do not move the auto-reader periodically or while it is being used. Connect the auto-reader to a secure and stable electrical mains outlet socket. Do not wet or heat the auto-reader. If liquid is spilled on the device, disconnect it and dry it immediately. Before turning it on, check that all reading positions are empty.

Using your auto-reader

- 1| Switch on your auto-reader by connecting the power supply's AC plug to the mains and then connect the plug at the other end of power supply to the rear of the auto-reader (See image 1).
- 2| Set the auto-reader local date and time running the Bionova® Wireless Assistant on a compatible mobile device.
- 3| Once the auto-reader turns on, the steps of the initialization sequence will start. Please wait until the initialization sequence has finished before using the auto-reader.
 - During the first step, the auto-reader will set and stabilize the incubation temperature. The temperature and stability indicator light will blink during this process. Once the working temperature is reached, the temperature and stability indicator will stop blinking.
 - In the final step, the auto-reader will stabilize the reading system for each incubation position. During this process, the Terragene® indicator light will blink. Once the step is complete, the Terragene® indicator light will stop blinking, and the initialization sequence has finished.

Note: The initialization sequence is carried out each time the auto-reader is turned on and may take between 10 and 30 minutes, depending on the current temperature of the auto-reader. Please, plan accordingly.

Note: Do not place an SCBI on an incubation position before the initialization sequence has finished.

- 4| Place an SCBI in an incubation position. The incubation will start automatically (See image 2)
- 5| Wait until the incubation has finished. The Position status indicator light will turn green to indicate a negative result, or red to indicate a positive result. If a positive result is found, press the Configuration & alarm cancellation button to silence the alarm. Do not remove the SCBI until a result has been produced. (See image 3)

- 6| Run Bionova® Cloud Web Application or the Bionova® Wireless Assistant Application to check the incubation result e-ticket.

Disposal: Dispose of the SCBIs in accordance with your country's waste-disposal guidelines.

Connectivity

To connect your auto-reader to a personal computer (PC) by USB, connect one end of the USB cable provided to the rear of the auto-reader and the other end to a USB port on your PC. Connect the auto-reader to a PC to use Bionova® Cloud Web Application. For more information, please refer to the User Manual. To access the extended features of your auto-reader using Wi-Fi or Bluetooth®, please download the Bionova® Wireless Assistant Application to your preferred mobile device and follow the steps detailed on the Application's tutorial.

Note: This auto-reader is only compatible with 2.4 GHz WiFi networks.

Note: For cybersecurity reasons, the auto-reader can only be connected to secure, password-protected Wi-Fi networks. It cannot be connected to open or public networks.

For a detailed overview on the connection process, please refer to the User Manual.

Connectivity credentials

During the connection process, the following credentials may be useful:

Device Internal Wi-Fi network credentials

Wi-Fi network name: PHOTON_[batch][serial number]

Wi-Fi network password: abcd1234

Bluetooth® credentials

Device Bluetooth® ID: PHOTON_[batch][serial number]

Auto-reader Incubation Programs

The following incubation program is available for your auto-reader.

Temperature selection	Time selection	Incubation program
60 °C	Instant (7 sec.)	Instant at 60 °C

Note: Please read and follow the SCBI's Instructions for Use before using it.

Troubleshooting

Refer to the User Manual for a detailed troubleshooting routine.

Please visit www.terragene.com/PHOTON or scan the QR code to access the auto-reader User Manual.

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