Discos con esporas Spore discs

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English

Composition

Each Bionova® paper disc is soaked with *Geobacillus stearothermophilus* ATCC 7953 spore population.

Product description

The Bionova® spore disc is specifically designed to control steam sterilization processes at 121 °C-134 °C. Spore discs should be incubated in Bionova® MC20 or MC1020 growth medium tube or in growth medium appropriate for *Geobacillus stearothermophilus* growth.

If sterilization process was not successful, MC20 or MC1020 growth medium will turn to yellow after incubation at 60±2 °C, thus indicating the existence of living *Geobacillus stearothermophilus* spores in the discs.

If sterilization process was correct, MC20 or MC1020 indicating medium will keep its original color. A last reading should therefore be done after 24 hours of growth medium incubation at 60±2 °C.

Warning!

Do not use spore discs to control radiation, ethylene oxide, dry heat or any sterilization processes other than steam.

Do not re-use the spore discs.

Storage

Store in a dark place and at temperatures between 10-30 °C, 30-80 % relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

Directions for use

1. Pack the Bionova® discs in appropriate package along with items for sterilization, according to recommended sterilization practices. Place the discs in those areas you a priori consider more inaccessible for sterilizing agent. A typical problematic area is the center of the load.

2. Sterilize as usual.

3. After sterilization process has finished, remove Bionova® discs from the package used for processing and incubation.

4. Transfer the spore discs with a sterile clamp, in sterility conditions to MC20 or MC1020 growth medium tube or to any other appropriate growth medium. **IMPORTANT!** Use latex gloves and chinstrap when transferring spore discs to growth medium tube.

WARNING! Avoid spore disc contact with any kind of surface including the outer part of growth medium tube.

5. Incubate spore disc at 60±2 °C.

IMPORTANT: Use a non-sterilized spore disc as a positive control each time a processed disc is incubated. The positive control ensure that correct incubation conditions were met.

6. Incubate processed spore discs along with the disc used as positive control during a maximum of 24 hours at 60±2 °C. It is advisable to make observations every 10 hours. Growth medium indicator color turning to yellow shows a failure on sterilization process. If after 24 hours no color change is visible on MC20 or MC1020 growth medium, the result is negative (i.e. the sterilization process was successful). The positive control indicator medium should show a yellow color change for the results to be valid. Record the positive ones and discharge them immediately as indicated below.

WARNING! Do not re-use sterilizer until spore discs growth result turns negative (MC20 or MC1020 growth medium containing processed spore disc remains in its original color).

Disposal

Dispose of the growth medium tubes already used to analyse Bionova® spore discs growth according to health regulations in your country. Positive MC20 or MC1020 growth medium tubes can be sterilized at 121 °C for at least 20 minutes, or at 132 °C for at least 15 minutes in gravity air displacement steam sterilizer, or at 134 °C for at least 10 minutes in vacuum assisted steam sterilizer.