
Chemical Tubes

Rev. 7 / 09.2021

CD210 - CD220 - CD250

Results reference guide

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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Explanation of Symbols

	Product designed for use in Steam sterilization cycles.		Batch number.
	Product designed for use in Dry Heat sterilization processes.		Manufacture Date.
	CE mark.		Expiration Date.
	Authorized representative in the European Community.		Manufacturer.

EN Chemical tubes

For sterilization processes control

Type 4 indicators according to ISO 11140-1:2014

Product description

Chemdye® CD210, CD220 and CD250 chemical tubes consist of an indicator solution contained in glass ampoules designed for controlling sterilization processes of liquid loads. Chemdye® CD210 and CD220 chemical tubes were

developed to monitor steam sterilization processes, while the chemical tube Chemdye® CD250 was designed to control dry heat sterilization processes.

•Initial color: **Red**

•End color: **Green**

Chemical tubes do not ensure proper sterilization of the product; they only make it possible to quickly determine if the critical parameters of the sterilization process for which they were designed have been reached.

The ampoules are sealed at both ends, ensuring the indicator solution contained does not come into contact with other elements during the sterilization cycle.

Control of steam sterilization processes

•CD210 (Black spot): Use in steam autoclaves at 121 °C for 15 minutes or at 126 °C for 10 minutes.

•CD220 (Yellow spot): Use in steam autoclaves at 134 °C autoclave for 3 to 3.5 minutes.

Control of dry heat sterilization processes

•CD250 (White spot): Use in sterilizing ovens at 160 °C for 120 minutes or 170 °C for 60 minutes or 180 °C for 35 minutes.

Instructions for use

1. Place the chemical tube in the center of the liquid load within the product to sterilize. For example: in the case of syringes, place the chemical tube inside the stem in order to ensure that the heat has penetrated even the most inner parts.

2. Sterilize as usual.

3. Once the sterilization process is complete, remove the chemical tube from the container and analyze the results by using the enclosed color guide. If the fluid remains red or amber, materials must be re-sterilized using a new tube.

IMPORTANT: Excessive exposure to high temperatures/long times may cause more subtle colors. Do not use at temperatures above 200 °C. If you use dry heat under 140 °C, the values required to assure sterility will not be achieved.

NOTE: It is not recommended that CD250 be exposed for periods exceeding double the suggested amount. After this amount of time the color of the solution returns to indicate that the sterilization cycle was not carried out correctly. This response serves as warning to the user, ensuring that materials are not exposed unnecessarily to stress caused by high temperatures, decreasing their useful life.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Precautions

Inspect the integrity of the chemical tube before use.

After the sterilization process is completed, always wear protective glasses and use gloves or tongs when handling the chemical tubes.

Allow the product to cool before handling.

Do not shake.

Do not use to control other sterilization processes than those indicated for each product.

Do not break the glass tube.

Do not reuse.

Storage

Keep away from light and refrigerate at a temperature between 2-15 °C, HR 30-80 %. Avoid freezing.

Shelf life

Chemical tubes have an expiration date of 2 years from the date of manufacture when stored at recommended storage conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Transport

Transport at a temperature below 30 °C for a maximum of 10 days in the case of CD210 and CD220 chemical tubes and for a maximum of 30 days in the case of CD250 chemical tubes.

Waste treatment

Processed tubes can be stored for recording the results, keeping them in a dark, cool place or disposing of them as sharps waste.

Color guide

Allows an instant visual check of the results. However, reference colors may not be exact.

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