CD83-ST

Rev. 0 / 06.2021

Chemical Indicators

Chemical Indicator
Indicador Químico

Failures
Fallas

Processed and correct
Procesado y correcto

EN

Chemical Indicator

For monitoring disinfection processes of surfaces by means of UV-C radiation

It has been demonstrated that many microorganisms can be inactivated by continuous UV-C radiation. UV-C radiation is absorbed by DNA and RNA, causing changes in their structures and preventing their replication. The radiation dose required to achieve a microbial load reduction depends not only on the particular microorganism, but also on additional factors like radiation intensity (power) of the emitting source, distance between emitting source and surface, and time of exposure. Furthermore, surfaces that only receive radiation by reflection, get lower doses that may not be sufficient to inactivate the pathogenic microorganisms present. Considering all this, it is of utmost importance to control that all surfaces to be disinfected, receive the minimum radiation dose required. The Chemdye® CD83-ST chemical indicator has been specially developed

to meet this need, having been calibrated to change its color due to exposure to UV-C light at doses greater than 50 mJ/cm².

Product description

The chemical indicator Chemdye® CD83-ST has been developed especially to monitor the surface disinfection processes that use UV-C as disinfectant agent. It consists of a self-adhesive polypropylene label containing reactive ink and a color reference for the final color change. The reactive indicating ink has been calibrated to experience different color changes depending on the accumulated radiation dose received. **IMPORTANT:** the reference color is illustrative and corresponds to a radiation of 50 mJ/cm².

The technology of the label enables its traceability by QR code, which allows its use in the DCS disinfection certification system.

The traceability QR code for the Disinfection Certificate provides information about the status of the label, indicating whether it has already been used in a Disinfection Process, the results, the data of the disinfection operator, the device used and the place, date and time in which it was subject to the procedure.

IMPORTANT: DCS CD83-ST indicators ensure that a dose of UV-C radiation has been delivered in an environment or surface, if its color changes to the reference color or more reddish. These indicators were designed for disinfection and do not ensure the cleanliness of the surface, which must be prior to the disinfection procedure.

Characteristics

29 x 22.5 mm self-adhesive polypropylene label, printed with indicator ink, for use in the Chemdye® CG3-N manual labeller (see instructions for use attached to the product).

Indicators 100 % free from toxic heavy metals

Smart Disinfection Certification System (DCS)

The CD83-ST labels have been designed as part of a brand new Certification Disinfection System. Thanks to its exclusive design, the DCS (Disinfection Certification System), enables the systematization and recording of the disinfection process by using the traceability given by its design which includes an unique QR code CD83-ST label.

Through its DCS App, the DCS uses a criterion based on artificial intelligence (IA) to guide the disinfection system that includes both good practices, suggestions and standardized disinfection protocols and makes it possible to interpret the result of the exposure to the disinfectant agent through direct reading with a Smartphone. For further information, please visit www.terragen.com/dcs-mobile-app/.

Instructions for use

1. Place the Chemdye® CD83-ST indicator on the surface where the disinfection level is to be checked according to the DCS protocol available in the DCS app.

NOTE: Recommended use: 4 Chemdye® CD83-ST chemical indicators per m². The indicators should be placed not higher than 1.5 m above ground level.

2. Run a normal disinfection cycle.

IMPORTANT: the label was design to be used only with UV-C radiation. The use of other disinfection methods and/or exposure to natural or artificial light may damage it.

3. Once the disinfection process is completed, compare the final color of each Chemdye® CD83-ST chemical indicator against the reference color strips printed in the same label or use the DCS App for the reading process. The final color it is related to the accumulated radiation dose received. If the final color is similar or darker than the reference, it means that the dose for which the label was designed has been reached (see the reference color chart).

NOTE: Chemdye® CD83-ST chemical indicators, both those exposed and those that have not been used, should be kept protected from light. Unexposed Chemdye® CD83-ST Chemical Indicators should be removed from their original packaging only just prior to running the disinfection

process and kept away from light and moisture.

Precautions

Do not expose to both sunlight and/or natural light as it may cause alterations in the reactivity of the indicator depending on the exposure time.

Exclusive use with UV-C light. Use with disinfectants or other sources of radiation, such as direct sunlight, could give erroneous results.

Do not reuse chemical indicators.

Product storage and expiration date

Store the Chemdye® CD83-ST in the bag in a dry place, away from light, at a temperature between 10 °C and 30 °C, at a relative humidity between 30 % and 80 %.

Do not wet.

Do not store near disinfectant/sterilizing agents.

Chemical indicators have an expiration date of 5 years from the date of manufacture if stored under recommended conditions.

Do not use the Chemdye® CD83-ST chemical indicator after its expiration date.

If the indicators are not used over a long time, the indicator should be removed from the automatic labeler and placed back into the original bag.

Waste treatment

Discard used chemical indicators according to your country's healthcare and safety regulations.

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Terragen S.A.