Helix Process Challenge Device

CDPCD2X025/P | CDPCD2X050/P | CDPCD2X075/P |
CDPCD2X10/P | CDPCD2X12/P | CDPCD2X15/P | CDPCD2X30/P |
CDPCD2X45/P | CDPCD2X025/M | CDPCD2X050/M |
CDPCD2X075/M | CDPCD2X10/M | CDPCD2X15/M |
CDPCD2X30/M | CDPCD2X45/M

EN

Helix Process Challenge Device

Product Description

Chemdye® Helix-PCD devices have been designed to evaluate the penetration capability of the sterilizing agent (Steam, Ethylene Oxide, Formaldehyde and Hydrogen Peroxide). For a successful sterilization, the sterilizing agent must contact all external and internal surfaces of the item to be sterilized. To achieve this, air and other non-condensable gases (NCG) must be evacuated from the sterilizer chamber prior to the sterilant injection. Chemdye® Helix-PCD devices simulate hollow instruments of different levels of complexity allowing testing of the penetration capability of the sterilizing agent.

Characteristics

Helix-PCD Chemdye® devices consists of a metal or plastic capsule of great resistance and durability, which is able to hold a chemical indicator designed for this system. The capsule has a screw cap which allows the introduction of the indicator and is connected to a high resistance plastic tube of variable length (for options available see our catalog) of 2 mm internal diameter. The plastic tube is helically arranged and it must not be released when using the product. A connector coated and reinforced by an external connector, holds the tube and capsule joint. Chemdye® Helix-PCD devices with the appropriate indicator must be processed through within the load. After the sterilization cycle, the result must be observed. Homogeneous color change of the chemical indicator ink means the sterilization process was successful and the load can be securely released.

Handling information

- 1- Choose an appropriate chemical indicator for the sterilization cycle to evaluate.
- 2- Fold the indicator strip in the middle through the punched line, the reactive ink must be facing the inside. Avoid folding the printed zone.
- 3- Place the chemical in the screw cap pin (metal capsule) or in the screw cap cavity (plastic capsule). Put the folded end first.
IMPORTANT: Avoid contact between chemical indicator and capsule inner surface, since the results can be altered.
- 4- Put the capsule cap by applying moderate pressure.
- 5- Place the system inside the fabric bag and then onto a stainless steel tray. Put the tray in the sterilization chamber near its door.
IMPORTANT: Do not use the fabric bag for Hydrogen Peroxide, Formaldehyde and Ethylene Oxide processes.
- 6- Run the desired sterilization cycle. For Ethylene Oxide processes, cycles shorter than 90 min or with less than 300 mg/l should not be evaluated with these devices.
- 7- Once the sterilization process is finished, remove the Chemdye® Helix-PCD from the sterilization chamber. **CAUTION:** For Steam processes, gloves should be worn. Also, allow 15 min for the load to cool down in order to avoid burns. For EO, Formaldehyde or Hydrogen Peroxide processes, gloves, glasses and mask should be worn to avoid contact with the agent remains. For EO, a 24-hour ventilation period or a time period suggested by the manufacturer, should be taken into consideration.
- 8- Remove the indicator and evaluate the result.
- 9- Record results by adhering the self-adhesive strip on your register book.

Chemical Indicator

Choose the appropriate chemical indicator according to the sterilization process to be evaluated.

Cycle	Indicator
Steam (3.5 min 134 °C)	PCD-A-3.5Y
Steam (5.3 min 134 °C 15 min. 121 °C)	PCD-A-5.3Y
Steam (7.0 min 134 °C 20 min. 121 °C)	PCD-A-7.0Y
Steam (9.0 min 134 °C)	PCD-A-9.0Y
Steam (18 min 134 °C)	PCD-A-18.0Y
Bowie-Dick (3.5 min. 134 °C)	PCD-A-3.5BD
Ethylene Oxide	PCD-A-E1
Formaldehyde	PCD-A-F1
Hydrogen Peroxide	PCD-A-P1

Interpretation of results

If the 4 bars have changed to the reference color, the sterilization cycle was successful. If one or more bars did not reach the final color, the load cannot be considered sterile.

Additional technical information

In each cycle, the user should check that the connector and the plastic tube are in perfect conditions. **NOTE:** In case the connector is damaged and replacement is needed, contact a sales representative and request a replacement. In case the tube is blocked or has leakages, the complete device should be replaced.

CDPCD2X15/P and CDPCD2X15/M (internal diameter 2 mm, tube length 1500 mm) devices are designed according to EN 867-5 standard. They can be used for performance testing for small steam sterilizers type B as well as for performance testing of large sterilizers that are used for processing hollow loads according to EN 285.

Chemdye® Helix-PCD devices and indicators are manufactured, calibrated and designed according to ISO 11140-1 standard.

Recommendations

Depending on its length, Chemdye® Helix-PCD devices show different degrees of challenge for a given sterilization cycle. Consequently, a non-homogenous change of the color of the indicator could be a consequence of choosing a device of inappropriate length. For example: for Plasma or Vaporized Hydrogen Peroxide sterilization cycles, it is recommended to use devices of the following lengths: 0.25, 0.5, 0.75 and 1.2 m. For routine load configurations, it is advised to evaluate the behavior of Chemdye® Helix-PCD devices of different degrees of challenge. After analyzing the results, for future cycles with the same load configuration, the most challenging device showing a complete color change in the 4 bars should be chosen.

Caution

Use Chemdye® Helix-PCD devices with metal capsules for Steam sterilization processes. Use Chemdye® Helix-PCD devices with yellow plastic capsules for EO, Formaldehyde and Hydrogen Peroxide sterilization processes; and Chemdye® Helix-PCD devices with grey plastic capsules for all processes listed above, plus Steam sterilization processes.
For EO processes it is not suggested to sterilize for less than 90 min and with gas concentration lower than 300 mg/l.
Terragene® does not guarantee proper functioning of the Chemdye® Helix-PCD device when used with indicators of other manufacturers.

Chemdye® Helix-PCD device results are only valid for the cycle run. Do not use them to draw conclusions about previous or future cycles.
Terragene® does not guarantee proper functioning of the Chemdye® Kit Helix-PCD when used for sterilization processes other than the process for which it has been designed.

Storage

Chemdye® Helix-PCD device and indicators should be stored in dry places, at temperatures between 10-30 °C, 30-80% RH, protected from light. Do not store close to sterilizing agents. Do not use indicators after their expiration date.

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