
## 3D Chemical Indicator

### Explanation of Symbols

Manufacture Date.  
Expiration Date.  
Manufacturer.

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### EN Chemical Indicators

For monitoring disinfection processes by continuous UV-C (254 nm) radiation

It has been demonstrated that many microorganisms can be inactivated by continuous UV-C radiation. UV-C radiation is absorbed by DNA and RNA, causing changes in their structures and preventing their replication. The radiation dose required to achieve a microbial load reduction depends not only on the particular microorganism, but also on additional factors like radiation intensity (power) of the emitting source, distance between emitting source and area to be disinfected, and time of exposure. Furthermore, shadowed areas which only receive radiation by reflection, get lower doses that may not be sufficient to inactivate the pathogenic microorganisms present. Considering all this, it is of utmost importance to control that all the areas of a room to be disinfected, receive the minimum radiation dose required.

### Product description

Chemdye® ChemSurf CDSUV-1 chemical indicator has been specially developed to control disinfection processes by continuous UV-C radiation. It consists of a pyramidal shape synthetic substrate of 10.5 cm side and 9.0 cm height printed with reactive indicator ink. The reactive indicator ink has been calibrated to experience different color changes depending on the accumulated radiation dose received, allowing to quickly and accurately control the efficacy of the disinfection process. Chemdye® ChemSurf CDSUV-1 chemical indicator exclusive design allows to dangle it in height and evaluate the accumulated radiation dose received in all directions. See forms of use in height and assemblies.

### Instructions for use

1. Mount as many three-dimension Chemdye® ChemSurf CDSUV-1 chemical indicators as areas to be verified, following the specified assembly instructions.

**NOTE 1:** it is recommended to put special emphasis on shadowed areas, which are those that do not receive radiation directly because something stands between them and the UV source, as they are particularly difficult to disinfect.

2. Leave the room and run a normal disinfection cycle.

3. Once the process has finished, see the final color of each side of each processed Chemdye® ChemSurf CDSUV-1 chemical indicator. If the final color matches the reference color, it means that the accumulated radiation dose received is equal to 200 mJ/cm<sup>2</sup>. If the final color is darker (redder) than the reference color, it means that the accumulated radiation dose received is higher than 200 mJ/cm<sup>2</sup>. If the final color is lighter (more orange) than the reference color, it means that the accumulated radiation dose received is lower than 200 mJ/cm<sup>2</sup>. The attached Documentation Chart allows to quickly register and analyze the obtained results. Chemdye® ChemSurf CDSUV-1 chemical indicator, has also a chart on one of its sides that allows to register its location in the room.

Confirm the color change with the reference color immediately or within 24 hours for best results.

**NOTE 2:** if Chemdye® ChemSurf CDSUV-1 chemical indicator is placed onto a surface, the side in contact with the surface will not experience any color change due to it will not have been exposed to the disinfectant agent. Chemdye® ChemSurf CDSUV-1 chemical indicator sides may experience different color changes as radiation received on every side might not be the same.

**NOTE 3:** exposed and non-exposed Chemdye® ChemSurf CDSUV-1 chemical indicators must be kept protected from light. Non-exposed Chemdye® ChemSurf CDSUV-1 chemical indicators must be removed from their original packaging ONLY right before running the disinfection process.

**NOTE 4:** if any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

### Precautions

Do not wet.

Do not store close to disinfectant or sterilizing agents.

### Storage and shelf life

Store in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80 %.

Chemdye® ChemSurf CDSUV-1 chemical indicator has an expiration date of 5 years from the date of manufacture if stored under recommended conditions. Do not use Chemdye® ChemSurf CDSUV-1 chemical indicator after its expiration date.

Once exposed, the indicator should be stored protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. The incorrect conservation of the chemical indicators may alter the color after exposure.

### Disposal

Discard used chemical indicators according to your country's healthcare and safety regulations.

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#### Figure 1

##### Instructions


1. Fold Chemdye® ChemSurf CDSUV-1 along the dotted lines, leaving the reactive indicating ink out.
2. Insert flap "A" in the slot corresponding to the same letter, as indicated in Fig.1. Repeat the same for all the flaps.
3. To dangle the indicator, pass the silicone seal through perforations of faces 2 and 3. Attach the self-adhesive hook to a suitable surface and dangle the pyramid from the seal.

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Figure 1 | Figura 1

**NOTE:** Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators. Images are reproduced for illustration purposes only, and do not necessarily match the article described in every detail.

Learn more about this product  
and check the digital color  
reference guide

Scan the QR code