Sequential integrator indicator

Quality Certification

STEAM

Product code

LOT

Temperature

Time (min.)

121 °C

128 °C

135 °C

*The curve allows to estimate the kill time within the temperature range informed.

According to ISO 11140-1:2014 standard

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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UNSAFE
SAFE
EXTENDED

EXTENDED condition

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

Explanation of Symbols

STEAM

Product designed for use with Steam sterilization cycles.

CE

CE mark.

EC REP

Authorized representative in the European Community.

LOT

Batch number.

Manufacture Date.

Expiration Date.

Manufacturer.

EN

Integrator Indicators

For Steam sterilization

Type 5 Indicator according to ISO 11140-1:2014

Product description
Integron® IT26 are three-level sequential integrating indicators designed to react inside loads exposed to Steam sterilization processes. The integration level consists of a pink ink developed to turn to black when a theoretical population of spores reaches its kill time, indicating that the **SAFE** condition is achieved.

Level 1 - **Blue: Exposure level**, it turns to black to indicate the exposure to the process but it does not assure that the conditions indicated in the integrator have been met (**UNSAFE** condition).
Level 2 - **Pink: Integration level**, it turns to black thus indicating that conditions indicated in the integrator have been met. These conditions are calibrated with the kill time of a 10⁶ theoretical *Geobacillus stearothermophilus* spore population, calculated in a BIER (Biological Indicator Evaluator Resistometer). This is called **SAFE** condition.
Level 3 - **Green: Over-integration level**, it turns to black when integration conditions indicated in the product have been overreached. This is called **EXTENDED** condition.

Characteristics
The indicator consists of a one-side-laminated paper strip of 60 mm x 25 mm.

Indications for use
Use integrator indicators for monitoring Steam sterilization processes. **SAFE** condition is reached in Steam sterilization processes as described in the quality certificate.

Precautions
Do not expose this product to EO, Dry Heat, or any sterilization processes other than Steam. Do not reuse the sterilizer until the integrator's color reaches or passes the **SAFE** condition.

Instructions of use
• Place an Integron® IT26 Integrator Indicator within each package, bag or tray that will be sterilized by Steam.
• In rigid containers place an Integron® IT26 Integrator Indicator in every corner, or at least two diagonally or opposite corners, to ensure that steam has reached the most inaccessible areas.
• After completing the process, analyze the results. If the **SAFE** integration condition has been reached or overreached, the sterilization process was successful.
NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Storage
Indicators should be stored in dry places, away from light, at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents.

Shelf life
Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date as well as the product batch number are specified on the package label and on each unit of product. Do not use indicators after their expiration date.
Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal
Discard along with cardboard waste according to your country's healthcare regulations.

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