PCD26-2

Rev. 7 / 09.2021

Process challenge device (PCD)

Quality certification

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STEAM

LOT

PCD26-2-RC Chemical Integrator

LOT

Stated Values (SV)

Temperature

Time

121 °C

min.

128 °C

min.

135 °C

min.

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

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Explanation of Symbols

STEAM

Product designed for use in Steam sterilization cycles.

CE

CE mark.

EC REP

Authorized representative in the European Community.

LOT

Batch number.

Manufacture Date.

Expiration Date.

Manufacturer.

Autoclave/Steam Cycles	
Model	Gravity Displacement
PCD26-2	121 °C
	30 minutes
	Dynamic Air Removal (Vacuum Assist)
	4 minutes
	132 °C

EN Process Challenge Device

For Steam sterilization

Indications for Use
United States
Terragene® Integron® PCD (PCD26-2) provides a defined challenge resistance against the claimed cycles shown below and demonstrated resistance equivalence to the AAMI/ANSI 16 towel pack. The device provides routine monitoring and sterilizer qualification testing steam sterilization processes. See Intended Use Table for reference.

Outside the United States
Integron® PCD26-2 Steam Process Challenge Device has been designed for a rapid and easy monitoring of dynamic air removal Steam sterilization processes at 132/135 °C ≥ 4 minutes and gravity air displacement processes at 121 °C ≥ 30 minutes.

Device Description
Integron® PCD26-2 Process Challenge Device (PCD) Test Pack has been designed to simulate a load to be sterilized and to pose a challenge to the sterilization process. It is used to evaluate the effective performance of the process by detecting inadequate air removal and steam penetration. Furthermore, it allows to perform the routine monitoring and periodic validation of the sterilizers (after repair, installation, relocation). Integron® PCD26-2 consists of a disposable pre-assembled package as outlined in ANSI/AAMI ST79 which contains in the middle a PCD26-2-RC self-adhesive Record Card printed with chemical integrator indicator (Type 5 according to ISO 11140-1:2014 standard) where sterilization cycle information can be recorded. The indicator ink was developed to change to a reference color when a theoretical spore population reaches its kill time, indicating integration conditions were reached. This condition is calibrated with the kill time of a 10⁶ *Geobacillus stearothermophilus* ATCC 7953 spore population, calculated in BIER (Biological Indicator Evaluator Resistometer). This pack also has a process indicator (Type 1 according to ISO 11140-1:2014 standard) printed on the outside of the package, which allows to quickly identify if it has been exposed to the sterilization process (color change from light blue to grey), and to distinguish between processed and unprocessed items.

Precautions
WARNING: Do not use PCD for monitoring Ethylene Oxide, Dry Heat, Formaldehyde or any sterilization process other than Steam.
WARNING: Place one or more PCDs in sterilizing hard-to-reach areas to ensure all areas of the chamber are sterilized. Evaluate all load configurations to ensure ALL hard-to-reach areas have been identified, and place a PCD in each of those locations.
WARNING: Do not use the sterilizer until the test result is correct.

Instructions for Use
1. Place the pack inside a normally loaded steam autoclave, in those areas which are considered most inaccessible for the sterilizing agent (e.g., the center of the load and areas near the door).
2. Run the sterilization cycle.
3. After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the test pack. **NOTE:** The color of the box may vary from the original after undergoing the sterilization cycle. This does

not represent a problem regarding the operation or quality of the product.
4. Check that the process indicator printed on box has changed color from light blue to grey. Open the test pack and remove the integrating indicator. **PRECAUTION:** Wear safety glasses and gloves when removing the indicator from the sterilized test pack.
5. Check the chemical integrator on the Record Card for correct exposure. Color change to the reference color confirms that the inside of the pack has been exposed to sufficient sterilization conditions. For chemical integrator reference color, please refer to Result Reference Guide. Otherwise, check the sterilization process.
6. Fill out the required information on the Record Card.
7. Discard the pack immediately.
NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Monitoring frequency
Follow facility Policies and Procedures which should specify a chemical indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As the best practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate indicator.

Storage
Store in a dark place, at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not freeze. Do not store near sterilizing agents or other chemical products.

Shelf Life
PCD26-2 has an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Do not use after expiration date. Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal
Discard the PCD after use along with paper waste, according to your country's healthcare and safety regulations.

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