
Results reference guide

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Chemical tubes

CD210 - CD220 - CD230 - CD240 - CD250

For sterilization processes control

Type 4 indicators according to ISO 11140-1:2014

Product description

Chemdye® CD210, CD220, CD230, CD240 and CD250 chemical tubes consist of an indicator solution contained in glass ampoules designed for controlling sterilization processes of liquid loads. Chemdye® CD210 and CD220 chemical tubes were developed to monitor steam sterilization processes, while Chemdye® CD230, CD240 and CD250 chemical tubes were designed to control dry heat sterilization processes.

- Initial color: **Red**
- End color: **Green**

Chemical tubes do not ensure proper sterilization of the product; they only make it possible to quickly determine if the critical parameters of the sterilization process for which they were designed have been reached.

The ampoules are sealed at both ends, ensuring the indicator solution contained does not come into contact with other elements during the sterilization cycle.

Control of steam sterilization processes

- CD210 (Black spot): Use in steam autoclaves at 121 °C for 15 minutes or at 126 °C for 10 minutes.
- CD220 (Yellow spot): Use in steam autoclaves at 134 °C autoclave for 3 to 3.5 minutes.

Control of dry heat sterilization processes

- CD230 (Green spot): Use in sterilizing ovens at 160 °C for 60 minutes or 170 °C for 30 minutes.
- CD240 (Blue spot): Use in convection sterilizing ovens at 180 °C for 12 minutes.
- CD250 (White spot): Use in sterilizing ovens at 160 °C for 120 minutes or 170 °C for 60 minutes or 180 °C for 35 minutes.

Instructions for use

-Place the chemical tube in the center of the liquid load within the product to sterilize. For example: in the case of syringes, place the chemical tube inside the stem in order to ensure that the heat has penetrated even the most inner parts. -Sterilize as usual.

-Once the sterilization process is complete, remove the chemical tube from the container and analyze the results by using the enclosed color guide. If the fluid remains red or amber, materials must be re-sterilized using a new tube.

IMPORTANT: Excessive exposure to high temperatures/long times may cause more subtle colors. Do not use at temperatures above 200 °C. If you use dry heat under 140 °C, the values required to assure sterility will not be achieved.

NOTE: It is not recommended that CD230, CD240 and CD250 be exposed for periods exceeding double the suggested amount. After this amount of time of the color of the solution returns to indicate that the sterilization cycle was not carried out correctly. This response serves as warning to the user, ensuring that materials are not exposed unnecessarily to stress caused by high temperatures, decreasing their useful life.

Precautions

Inspect the integrity of the chemical tube before use.

After the sterilization process is completed, always wear protective glasses and use gloves or tongs when handling the chemical tubes.

Allow the product to cool before handling.

Do not shake.

Do not use to control other sterilization processes than those indicated for each product.

Do not break the glass tube.

Do not reuse.

Storage and shelf life of the product

Keep away from light and refrigerate at a temperature between 2-15 °C, HR 30-80%. Avoid freezing.

Chemical tubes have an expiration date of two years from the date of manufacturing when stored under the recommended conditions. The expiration date and the product batch are specified on the package label. Do not use chemical tubes after their expiration date.

Transport

Transport at a temperature below 30 °C for a maximum of 10 days in the case of CD210 and CD220 chemical tubes and for a maximum of 30 days in the case of CD230/240/250 chemical tubes.

Waste treatment

Processed tubes can be stored for recording the results, keeping them in a dark, cool place or disposing of them as sharps waste.

Color guide

Allows an instant visual check of the results. However, reference colors may not be exact.

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