SKEO

Cintape®, Chemdye® and Integron®
Chemical Indicators + Bionova®
Biological Indicators + Bionova®
MiniBio Auto-Reader Incubator Kit

For Ethylene Oxide sterilization cycles control

Usage

Ethylene Oxide sterilization processes control.

Applicable Regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

- Bionova® MiniBio Auto-Reader Incubator: European Union Low Voltage Directive 2014/35/EU, Electromagnetic Compatibility Directive 2014/30/CE, RoHS Directive 2011/65/EU.
- Bionova® BT10, BT110, BT40 Biological Indicators and the culture medium Bionova® MC1030-2: ISO 11138-1:2017, ISO 11138-2:2017 and ISO 11138-4:2017; IRAM 37102-1:1999 and IRAM 37102-2:1999.
- Cintape® CT10, Chemdye® CD16, Integron® IT12 Chemical Indicators: ISO 11140-1:2014.

Classification

Class 1 medical devices, according to risk (Chemical and Biological Indicators).

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1 and PM 1614-4.

Characteristics

- Bionova® MiniBio Auto-Reader Incubator
 - 3 metal positions (38 mm deep, 9 mm diameter) to incubate Rapid/Super Rapid/Ultra Rapid Fluorescence Readout Self-contained biological indicators.
 - 1 ampoule crusher.
 - Special hole for external thermometer.
 - Incorporated top cap.
 - Dual Temperature System: Allows selecting between two different incubation temperatures (37 °C and 60 °C).
 - USB connection for registering results in a PC through Reading and Traceability Software.
 - Allows incubation of 3 indicators simultaneously, in different incubation times, at the same temperature.
 - Thermal Printer to record results.
 - Sound alarm for event indication.
 - Detection and automatic cancellation of the fluorescence reading of the biological indicator.
 - Fluorescence reading wavelength: ~ 460 nm.
 - Firmware update by the user.
 - Dimensions: 180 mm high, 154 mm wide and 105 mm deep.
 - Voltage range: 100 – 240 V AC.

Power: 28 W.

Frequency: 50-60 Hz.

- Biological Indicators

Bionova® BT10

Self-Contained Biological Indicator.

10^6 *Bacillus atrophaeus* ATCC 9372 spores per vial.

Incubation at (37 ± 2) °C for 48 hours (maximum).

Positive growth results (failure in the sterilization process): Culture medium turns to yellow.

Bionova® BT110

Self-Contained Biological Indicator.

$\geq 10^6$ *Bacillus atrophaeus* ATCC 9372 spores per vial.

Rapid Readout Fluorescence System. Final fluorescence reading is performed after 4-hour incubation at (37 ± 2) °C.

An optional visual pH color change confirmation could be made after incubation for 48 hours at (37 ± 2) °C.

If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at (37 ± 2) °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain blue after the incubation process.

7-day readout for visual confirmation by color change is optional and not intended to be routinely performed; it is an initial validation of the 4 hour-reading. Fluorescence results may be compared to the 7-day visual reading. NOTE: if 7-day readout is performed, a humidified environment will be required to avoid medium dry out.

Bionova® BT40

Spore Strip

10^6 *Bacillus atrophaeus* ATCC 9372 spores per strip.

The spore strip should be tested in a suitable culture medium at (37 ± 2) °C. We suggest the use of Bionova® MC1030-2 or MC1020-2 culture media to final reading in 48 hours.

Bionova® MC1030-2

Sterile growth medium is specially formulated for use with our Bionova® products (also suitable for other brand products).

Growth medium color: blue

Plastic tube: 68 mm high, 10 mm diameter.

Medium volume: 2 ml

For germination and growth of *Bacillus atrophaeus* ATCC 9372 and *Bacillus subtilis* ATCC35021 bacterial spores.

Incubation conditions: 48 hours (maximum) at (37 ± 2) °C.

Medium turns to yellow when spores grow.

- Chemical Indicators

Cintape® CT10

Type 1 Chemical Indicator.

Roll of self-adhesive tape printed with indicator ink..

Size: 50 m x 18 mm.

Initial color: orange.

Final color: green.

100 % Toxic Heavy Metal free

Chemdye® CD16

Type 4 Multivariable process Indicator.

Double punched card strips printed with indicator ink.

Size: 200 x 15 mm

Initial color: Púrpura/brown

Final color: Green

100 % Toxic Heavy Metal free

Integron® IT12

Type 5 Chemical Integrator.

Card printed with indicator ink, laminated on the back side.

Size: 60 mm x 25 mm

Level 1 – Red: Exposure level, it turns to yellow to indicate the exposure to the process but it does not assure that the conditions indicated in the integrator have been met (UNSAFE condition).

Level 2 – Purple/Brown: Integration level, it turns to green thus indicating that conditions indicated in the integrator have been met. These conditions are calibrated with the kill time of a 10^6 *Bacillus atrophaeus* spores population, calculated in BIER (Biological Indicator Evaluator Resistometer). This is called SAFE condition. Conditions: EO 600 mg/l, RH 60 %, 54 °C.

100 % Toxic Heavy Metal free

Environmental conditions during manufacture

Bionova® MiniBio Auto-Reader Incubator: T= 15-35 °C, RH= 35-85 %.

Bionova® Biological Indicators: T= 15-30 °C, RH= 30-80 %. Sterility conditions are necessary only during inoculation process performed in laminar flow.

Cintape®, Chemdye® and Integron® Chemical Indicators: T= 15-30 °C, RH= 30-80 %.

Storage conditions

T = 10-30 °C, RH= 30-80 %, keep in a dark place.

Transport conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent any risk for human health.

Shelf life

Bionova® MiniBio Auto-Reader Incubator: 5 years.

Bionova® Biological Indicators: 2 years.

Cintape®, Chemdye® and Integron® Chemical Indicators: 5 years.

Packing

1 unit of Bionova® MiniBio, 1 foil bag with 10 units of Bionova® BT10, 1 foil bag with 10 units of Bionova® BT110, 1 foil bag with 10 units of Bionova® BT40, 1 foil bag with 10 units of Bionova® MC1030-2, 1 foil bag with 1 roll of self-adhesive tape Cintape® CT10, 1 foil bag with 250 units of Chemdye® CD16, 1 foil bag with 250 units of Integron® IT12 and a catalogue per box.

Box information: product description, process for intended use, manufacturer information and data on box's label.

Labelling

On Bionova® MiniBio: same label on product and product's box: product code and description, batch number, serial number and manufacturer information.

On each Biological Indicators vial: 17.0 mm x 33.0 mm polypropylene label. Printed in black. 1.5 mm chemical indicator line, printed with Ethylene Oxide reactive ink (Color change to green). Product code, process for intended use, name of organism, batch number and expiration date.

On each Culture Medium vial: Label 17 mm x 33 mm printed in black. Product code, batch number, expiration date and bacterial strain that can be incubated.

On each foil bag: product code and description, process for intended use, presentation, applicable regulation, batch number, manufacture and expiration date, bacterial strain and population on Biological Indicators, color change on Chemical Indicators.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information

Bionova® MiniBio device has a button located on the incubator panel to select incubation temperatures at 37 ° C or 60 ° C. It also has a button for each of the positions of readout to set up a different time for each position. Bionova® Mini Bio incubator has an automatic readout system and a printer to record the results.

Read product's instructions for use thoroughly before use.

Precautions

Do not pour any liquid inside of Bionova® MiniBio incubator. Do not immerse into any liquid. Only for indoor use. Disconnect the power cord before cleaning. Do not use abrasive or corrosive cleaners or disinfectants. Make sure that the incubator is connected to a properly rated power cord.

Do not store the products near sterilizing agents.

Do not expose the products to Steam, Dry Heat, Radiation or any sterilization process other than Ethylene Oxide.

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