SKHP

**Kit Lectora-Incubadora Bionova®
MiniBio + Indicadores Biológicos
Bionova® + Indicadores Químicos
Chemdye® y Cintape®**

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For Plasma or Vaporized Hydrogen Peroxide sterilization cycles control

Usage

Plasma or Vaporized Hydrogen Peroxide sterilization processes control.

Applicable Regulation

- Bionova® MiniBio Auto-Reader Incubator: European Union Low Voltage Directive 2014/35/EU, Electromagnetic Compatibility Directive 2014/30/CE, RoHS Directive 2011/65/EU.
- Bionova® BT91, BT95, BT96 Biological Indicators: ISO 11138-1:2006; IRAM 37102-1:1999.
- Cintape® CT40, Chemdye® CD40, Chemdye® CD42, Chemdye® CD47 Chemical Indicators: ISO 11140-1:2014.

Classification

Class 1 medical devices, according to risk (Chemical and Biological Indicators).

Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1 and PM 1614-4.

Characteristics

- Bionova® MiniBio Auto-Reader Incubator
 - 3 metal positions (38 mm deep, 9 mm diameter) to incubate Rapid/Super Rapid/Ultra Rapid Fluorescence Readout Self-contained biological indicators.
 - 1 ampoule crusher.
 - Special hole for external thermometer.
 - Incorporated top cap.
 - Dual Temperature System: Allows selecting between two different incubation temperatures (37°C and 60°C).
 - USB connection for registering results in a PC through Reading and Traceability Software.
 - Allows incubation of 3 indicators simultaneously, in different incubation times, at the same temperature.
 - Thermal Printer to record results.
 - Sound alarm for event indication.
 - Detection and automatic cancellation of the fluorescence reading of the biological indicator.
 - Fluorescence reading wavelength: ~ 460 nm.
 - Firmware update by the user.
 - Dimensions: 180 mm high, 154 mm wide and 105 mm deep.
 - Voltage range: 100 – 240 V AC.
 - Power: 28 W.
 - Frequency: 50-60 Hz.

- Bionova® Biological Indicators

BT91

10⁵ or 10⁶ *Geobacillus stearothermophilus* ATCC 7953 spores per vial.
Positive growth results (failure in the sterilization process): culture medium turns to yellow, incubation at 60 °C for 24 hours (maximum).

BT95

Rapid Readout Fluorescence System.

Final fluorescence reading is performed after 2-hour incubation at 60 °C.

BT96

Super Rapid Readout Fluorescence System.

Final fluorescence reading is performed after 30-minute incubation at 60 °C.

≥ 10⁶ *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

An optional visual pH color change confirmation could be performed after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed; it is an initial validation of the 30-minute reading. Fluorescence results may be compared to the 7-day visual reading.

NOTE: If 7-day readout is performed, a humidified environment will be required to prevent medium from drying out.

D-Value at 2 mg/l H₂O₂, 50 °C.

- Chemical Indicators

Cintape® CT40

Type 1 Chemical Indicator.

Roll of self-adhesive tape printed with indicator ink.

Size: 50 m x 19 mm.

Initial color: purple

Final color: green

100 % Toxic Heavy Metals free

Tape surface treated with covering materials that allow an easy and continuous writing. Adhesive capable of adhering to wrap cloth, paper and plastic. Designed with advanced technologies to ensure that the tape remains adhered to the surface after sterilization and leaves no trace of adhesive on the surface of the materials after removal.

Chemdye® CD40
Type 4 Chemical Indicator.
Strips printed with indicator ink.
Size: 105 x 18 mm
Initial color: Púrples
Final color: Green
100 % Toxic Heavy Metals free

Chemdye® CD42
Type 1 Chemical Indicator.
Varnished strips printed with indicator ink.
Size: 105 x 18 mm
Initial color: Red
Final color: Yellow
100 % Toxic Heavy Metals free.

Chemdye® CD47
Type 1 Chemical Indicator.
Self-adhesive label (dot) printed with indicator ink.
Size: 15 mm of diameter
Initial color: purple
Final color: green
100 % Toxic Heavy Metals free

Environmental conditions during manufacture

Bionova® MiniBio Auto-Reader Incubator: Room temperature 15-35 °C, RH 35-85 %.

Bionova® Biological Indicators: Room temperature 15-30 °C, RH 30-80 %.

Sterility conditions are necessary only during inoculation process performed in laminar flow.

Chemical Indicators: 15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place.

Transportation conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf-life

Bionova® MiniBio Auto-Reader Incubator: 10 years.

Bionova® Self-Contained Biological Indicators: 2 years.

Chemical Indicators: 5 years.

Packing

1 unit of Bionova® MiniBio, 1 foil bag with 10 units of Bionova® BT91, 1 foil bag with 10 units of Bionova® BT95, 1 foil bag with 10 units of Bionova® BT96, 1 unit of Cintape® CT40, 1 foil bag with 250 units of Chemdye® CD40, 1 foil bag with 250 units of Chemdye® CD42, 1 foil bag with a roll of 1000 labels of Chemdye® CD47 and a catalogue per box.

Box information: product description, process for intended use, manufacturer information and data on box's label.

Packaging information on each product: product description, manufacturer information and data on box's label.

Labelling

On Bionova® MiniBio: same label on product and product's box: product code and description, batch number, serial number and manufacturer information.

On each Biological Indicators vial: 17.0 mm x 33.0 mm polypropylene label. Printed in black. 1.5 mm chemical indicator line, printed with Hydrogen Peroxide reactive ink (Color change to green). Product code and batch number, process for intended use and bacterial strain.

On each foil bag: product code and description, process for intended use, product presentation, regulation applied to the product, manufacturer information, batch number, manufacture and expiration date, bacterial strain and population on Biological Indicators, operating instructions and classification according to regulation on Chemical Indicators.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information (optional)

Bionova® MiniBio device has a button located on the incubator panel to select incubation temperatures at 37 ° C or 60 ° C. It also has a button for each of the positions of readout to set up a different time for each position.

Bionova® Mini Bio incubator has an automatic readout system and a printer to record the results.

Read product's directions for use thoroughly before use.

NOTE: When necessary, limits of each product's parameters should be specified.

Description of identified and/or predictable risks and situations which could represent a risk for the product

Do not pour any liquid inside of MiniBio incubator.

Do not immerse into any liquid.

Only for indoor use.

Disconnect the power cord before cleaning.

Do not use abrasive or corrosive cleaners or disinfectants.

Make sure that the incubator is connected to a properly rated power cord.

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Hydrogen Peroxide.

Once the Biological Indicators are used it is recommended to sterilize at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilizer, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.