SKST

**Kit Lectora-Incubadora Bionova®
MiniBio + Indicadores Biológicos
Bionova® + Indicadores Químicos
Chemdye®, Cintape® e Integron®**

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SKST

Chemdye®, Cintape®, Integron®
Chemical Indicators + Bionova®
Biological Indicators + Bionova®
MiniBio Auto-Reader Incubator Kit

For Steam sterilization cycles control

Usage

Steam sterilization processes control.

Applicable Regulation

- Bionova® MiniBio Auto-Reader Incubator: European Union Low Voltage Directive 2014/35/EU, Electromagnetic Compatibility Directive 2014/30/CE, RoHS Directive 2011/65/EU.
- Bionova® BT20, BT220, BT222 and BT224 Biological Indicators: ISO 11138-1:2006 and ISO 11138-3:2006; IRAM 37102-1:1999 and IRAM 37102-3:1999.
- Cintape® CT22, Chemdye® CD29, Integron® IT26-1YS, Integron® IT27-3YS and Integron® IT26-C Chemical Indicators: ISO 11140-1:2014.
- Chemdye® BD125X/2 Chemical Indicators: ISO 11140-1:2014 and ISO 11140-4:2007.

Classification

Class 1 medical devices, according to risk (Chemical and Biological Indicators).

Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1 and PM 1614-4.

Characteristics

- Bionova® MiniBio Auto-Reader Incubator
 - 3 metal positions (38 mm deep, 9 mm diameter) to incubate Rapid/Super Rapid/Ultra Rapid Fluorescence Readout Self-contained biological indicators.
 - 1 ampoule crusher.
 - Special hole for external thermometer.
 - Incorporated top cap.
 - Dual Temperature System: Allows selecting between two different incubation temperatures (37°C and 60°C).
 - USB connection for registering results in a PC through Reading and Traceability Software.
 - Allows incubation of 3 indicators simultaneously, in different incubation times, at the same temperature.
 - Thermal Printer to record results.
 - Sound alarm for event indication.
 - Detection and automatic cancellation of the fluorescence reading of the biological indicator.
 - Fluorescence reading wavelength: ~ 460 nm.
 - Firmware update by the user.
 - Dimensions: 180 mm high, 154 mm wide and 105 mm deep.
 - Voltage range: 100 – 240 V AC.
 - Power: 28 W.
 - Frequency: 50-60 Hz.

- Bionova® Biological Indicators

BT20

10^5 or 10^6 *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

Positive growth results (failure in the sterilization process): culture medium turns to yellow, incubation at 60 °C for 24 hours (maximum).

BT220

Rapid Readout Fluorescence System.

Final fluorescence reading is performed after 3-hour incubation at 60 °C.

BT222

Super Rapid Readout Fluorescence System.

Final fluorescence reading is performed after 1-hour incubation at 60 °C.

BT224

Ultra Rapid Readout Fluorescence System.

Final fluorescence reading is performed after 20-minute incubation at 60 °C.

$\geq 10^6$ *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

An optional visual pH color change confirmation could be made after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process was successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed.

Fluorescence results may be compared to the 7-day visual reading.

NOTE: If 7-day readout is performed, a humidified environment will be required to avoid medium to dry out.

D-Value: Not lower than 1.5 minutes at 121 °C. Another D-value is informed at 132°C and/or 135°C.

Z-Value: Not lower than 6 °C.

- Chemical Indicators

Cintape® CT22

Type 1 Chemical Indicator.

Roll of Super Adherent self-adhesive tape printed with indicator ink.

Size: 50 m x 19 mm.

Initial color: yellow

Final color: dark brown/black

Tape surface treated with covering materials that allow an easy and continuous writing. Super adhesive capable of adhering to wrap cloth, paper and plastic. Designed with advanced technologies to ensure that the tape remains strongly adhered to the surface after sterilization and leaves no trace of adhesive on the surface of the materials after removal.

Chemdye® CD29

Type 4 Multi-variable Chemical Indicator.

One-side laminated punched double card strips printed with indicator ink.

Size: 200 x 15 mm.

Initial color: yellow.

Final color: black.

100 % Toxic Heavy Metals free.

Integron® IT26-C

Type 5 Integrator Indicator.

Integrator Indicator consists of a paper wick and a steam and temperature sensitive chemical pellet contained in a paper/film/foil laminate. The chemical pellet melts and migrates as a dark bar along the paper wick. The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met. The extent of migration depends on steam quality, time, and temperature. The ACCEPT result is reached when a theoretical spore population reaches its kill time, indicating integration condition has been reached. This condition is calibrated with the kill time of a 10^6 *Geobacillus stearothermophilus* spores population, calculated in BIER (Biological Indicator Evaluator Resistometer).

Size: 52 mm x 21 mm.

Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

Integron® IT26-1YS

Type 5 Integrator Indicator.

One-side-laminated strips printed with indicator ink.

Size: 72 mm x 23 mm.

Initial color: yellow

Final color: black

100 % Toxic Heavy Metals free

Integration condition calibrated with the kill time of a 10^6 *Geobacillus stearothermophilus* spores population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

The product ensures an adequate control of the efficacy of sterilization processes (temperature, time, quality of steam).

Integron® IT27-3YS

Type 6 Emulator Indicator.

One-side-laminated strips printed with indicator ink.

Size: 72 mm x 23 mm.

Initial color: yellow

Final color: black

100 % Toxic Heavy Metals free

The product ensures an adequate control of the efficacy of sterilization processes (temperature, time, quality of steam).

Chemdye® BD125X/2

Type 2 Chemical Indicator

Initial Color: yellow

Final Color: dark brown/black (should change homogeneously in the whole sheet)

Size: 16 mm x 125 mm x 125 mm

Assembly: 19 permeable barriers - 1 Glassine paper - 1 B&D Test Sheet - 19 permeable barriers. Wrapped with crepe paper.

Permeable barriers: 200 g/m².

B&D Test Sheet: 240 g/m².

100 % Lead Free

Air removal or steam quality problems will be shown as heterogeneous color change into dark brown/black on the test sheet, such as lighter areas in the center or a non-uniform color change.

Environmental conditions during manufacture

Bionova® MiniBio Auto-Reader Incubator: Room temperature 15-35 °C, RH 35-85 %.

Bionova® Biological Indicators: Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during inoculation process performed in laminar flow.

Chemical Indicators: 15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place.

Transportation conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf-life

Bionova® MiniBio Auto-Reader Incubator: 10 years.

Bionova® Self-Contained Biological Indicators: 2 years.

Chemical Indicators: 5 years.

Packing

1 unit of Bionova® MiniBio, 1 foil bag with 10 units of Bionova® BT20, 1 foil bag with 10 units of Bionova® BT220, 1 foil bag with 10 units of Bionova® BT222, 1 foil bag with 10 units of Bionova® BT224, 1 unit of Cintape® CT22, 1 foil bag with 250 units of Chemdye® CD29, 1 foil bag with 100 units of Integron® IT26-C, 1 foil bag with 200 units of Integron® IT26-1YS, 1 foil bag with 200 units of Integron® IT27-3YS, 2 units of Chemdye® BD125X/2 and a catalogue per box.

Box information: product description, process for intended use, manufacturer information and data on box's label.

Packaging information on each product: product description, manufacturer information and data on box's label.

Labelling

On Bionova® MiniBio: same label on product and product's box: product code and description, batch number, serial number and manufacturer information.

On each Biological Indicators vial: 17.0 mm x 33.0 mm polypropylene label. Printed in black. 1.5 mm chemical indicator line, printed with Steam reactive ink (Color change to brown). Product code and batch number, process for intended use and bacterial strain.

On each foil bag: product code and description, process for intended use, product presentation, regulation applied to the product, manufacturer information, batch number, manufacture and expiration date, bacterial strain and population on Biological Indicators, operating instructions, classification according to regulation and stated values (SV) on Chemical Indicators.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information (optional)

Bionova® MiniBio device has a button located on the incubator panel to select incubation temperatures at 37 ° C or 60 ° C. It also has a button for each of the positions of readout to set up a different time for each position.

Bionova® Mini Bio incubator has an automatic readout system and a printer to record the results.

Read product's directions for use thoroughly before use.

NOTE: When necessary, limits of each product's parameters should be specified.

Description of identified and/or predictable risks and situations which could represent a risk for the product

Do not pour any liquid inside of MiniBio incubator.

Do not immerse into any liquid.

Only for indoor use.

Disconnect the power cord before cleaning.

Do not use abrasive or corrosive cleaners or disinfectants.

Make sure that the incubator is connected to a properly rated power cord.

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

Once the Biological Indicators are used it is recommended to sterilize at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilizer, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.