CDPCD2X30/M

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CDPCD2X30/M

Helix Process Challenge Device

To evaluate air removal and penetration of sterilizing agent in hollow loads.

Usage

To evaluate air removal and penetration of steam in hollow loads.

Conditions:

STEAM: 3.5 min., 134 °C // 5.3 min., 134 °C // 7.0 min., 134 °C // 9.0 min., 134 °C // 18.0 min., 134 °C.

Applicable regulation

More demanding requirements than EN 867-5.

Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012.

Characteristics

Metal capsule of great resistance and durability, to hold a chemical indicator designed for this system. The capsule is connected to a high resistance plastic tube of 3.0 m length and 2 mm internal diameter. Two special silicone connectors coated by external connector.

Chemical indicators designed for Helix-PCD system: PCD-A-3.5 / PCD-A-5.3 / PCD-A-7.0 / PCD-A-9.0 / PCD-A-18 (STEAM)

NOTE: If 4 bars change to the final color, the result is successful, the sterilizer works properly. If one or more bars have not turned to the final color, this shows a failure in sterilizing agent penetration, therefore an equipment malfunction.

Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep out of direct light.

Transportation conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

Metal capsule: 500 cycles. Plastic helix: 250 cycles (replaceable at very low cost).

Packing

1 Helix-PCD + 1 Cotton bag

Packing information: product code and description, manufacturer information and data on label.

Labelling

Only on product's box: product code and batch number.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information (optional)

Read product's directions for use thoroughly before use.

Note: When necessary, limits of each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to chemical agents.