Infection Control Division | Chemical Indicators

CT22 Indicator Tape

For steam sterilization processes

Usage

United States

Chemical process indicator tape intended for monitoring the efficacy of steam sterilization processes. Intended use cycles:

Gravity Displacement: 121 °C for 30 minutes, 132 °C for 15 minutes, 132 °C for 25 minutes, 135 °C for 10 minutes.

Dynamic Air Removal (Vacuum Assist): 132 °C for 4 minutes, 135 °C for 3 minutes.

Outside the United States

Sterilization package sealing tape with process indicator for monitoring steam sterilization processes. Designed for gravity displacement and dynamic air removal (vacuum assisted) sterilizers operating from 121 °C to 135 °C.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN, ISO 13485:2016.

ISO 11140-1:2014.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Classification

Class 1, according to risk (ANMAT).

FDA 510(k)

K191021.

Characteristics

Type 1 Chemical Indicator.

Roll of self-adhesive tape printed with indicator ink.

Size: 50 m x 18 mm.

Initial color: **yellow**.

Final color: **dark brown/black**.

100 % Toxic Heavy Metal free.

Tape surface treated with covering materials that allow an easy and continuous writing. Capable of adhering to wrap cloth, paper and plastic. Designed with advanced technologies to ensure that the tape remains adhered to the surface after sterilization and leaves no trace of adhesive on the surface of the materials after removal.

Infection Control Division | **Chemical Indicators**

CT22 Indicator Tape

For steam sterilization processes

Environmental conditions during manufacture

T = 15-30 °C, RH = 30-80 %.

Storage conditions

T = 10-30 °C, RH = 30-80 %, keep out of direct light.

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 10 - 30 °C.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

Acceptable excursions during transport are as follows:

T = 2-10 °C and T = 30-37 °C for no longer than seven days. No restrictions for relative humidity.

Shelf life

5 years.

Packing

1 roll of 50 m x 18 mm.

Packing information: product description, regulation, storage conditions, presentation, manufacturer information and data on packing's label.

Labeling

On roll's cone: process for intended used, classification according to regulation, batch number, expiration date and manufacturer.

On product's packing: product code and description, process for intended used, classification according to regulation, color change, tape dimensions, batch number, expiration date, bar code and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other relevant information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.

Do not expose this product to ethylene oxide, dry heat or any sterilization processes other than steam.

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