
IT26-AD

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
IT26-AD

Three level Integrator Indicator

For Steam sterilization processes.

Usage

Control of Steam sterilization processes.

Applicable Regulation

ISO 11140-1:2014

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management Systems standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-4.

Characteristics

Type 5 Integrator Indicator.

Self-adhesive one-side-laminated card printed with indicator ink.

Size: 60 mm x 25 mm.

Level 1 – Blue: Exposure level, it turns to black to indicate exposure to the process but it does not assure that the conditions indicated in the integrator have been met (UNSAFE condition).

Level 2 – Pink: Integration level, it turns to black thus indicating that the conditions indicated in the integrator have been met. These conditions are calibrated against the kill time of a 10^6 *Geobacillus stearothermophilus* spore population, calculated in a BIER (Biological Indicator Evaluator Resistometer). This is called SAFE condition. Conditions: saturated steam at 121 °C, 128 °C, 135°C.

Level 3 – Green: Over-integration level, it turns to black when integration conditions indicated in the product have been exceeded. This is called EXTENDED condition.

Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30°C, RH 30-80 %, keep out of direct light.

Transportation conditions

Storage conditions should be strictly followed. Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

5 years

Packing

250 strips per aluminium foil bag.

Packing information: product code and description, process for intended use, operating instructions, product presentation, classification according to regulation, batch number, manufacturer information and expiration date.

Weight per bag: 106, 3 g

Labelling

Information only on product's package.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information (optional)

Read product's directions for use thoroughly before use.

Note: When necessary, limits for each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat or any sterilization processes other than Steam.

Do not reuse the sterilizer until the indicator has changed into the color indicated in the product.