Infection Control Division | **Chemical Indicators**

KPCD222-C Steam Test Pack PCD Kit

Process Challenge Device for Steam sterilization processes

Usage

United States

Terragene® Bionova® PCD222-C provides a defined challenge resistance against the claim cycles shown below and also demonstrated resistance equivalence to the ANSI/AAMI 16 towel pack. The device provides routine monitoring and sterilizer qualification testing steam sterilization processes.

Autoclave/Steam Cycles		
Gravity Displacement		
121 °C	132 °C	135 °C
30 min	15 min , 25 min	10 min
Dynamic Air Removal (Vacuum Assist)		
132 °C		135 °C
4 min		3 min
Fluorescence Read Time and pH Color Change		
Fluorescence Read Time		pH Color Change
1 h		48 hs

Outside the United States

Bionova® PCD222-C Steam Process Challenge Device has been designed for quick and easy monitoring of dynamic air removal and gravity air displacement Steam sterilization processes at 132/135 °C ≥ 4 minutes and at 121 °C ≥ 30 minutes.

Applicable regulation

Designed under Quality Management Systems standards ISO 13485:2016/NS-EN, ISO 13485:2016.

ISO 11140-1:2014, ISO 11138-1:2017, ISO 11138-3:2017, ANSI/AAMI ST79.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Classification

Class 1, according to risk (ANMAT).

FDA 510(k)

K163646.

Infection Control Division | **Chemical Indicators**

KPCD222-C Steam Test Pack PCD Kit

Process Challenge Device for Steam sterilization processes

Characteristics

The kit consists of:

- PCD222-C: Porous cards system holding:

- A Self-contained Biological Indicator Bionova® BT222.

- A Type 5 chemical indicator Integron® IT26-C (moving-front integrator).

- A self-adhesive Record Card, PCDBI-C-RC, where data on sterilization cycle may be written.

- Self-contained Biological Indicator Bionova® BT222 (to be used as a positive control).

Process Indicator:

Initial color: **blue**.

Final color: **dark grey or black**.

IT26-C Integrator Indicator (within the challenge device):

The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met.

Integration condition is calibrated against the kill time of a 10^6 *Geobacillus stearothermophilus* ATCC® 7953 spore population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

BT222 Self-contained Super Rapid Biological Indicator (within the challenge device):

$\geq 10^6$ *Geobacillus stearothermophilus* ATCC® 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

Final results: by fluorescence after incubation for 1 hour at 60 °C (sensitivity 97%). A visual confirmation of the result is optional (color change of the culture medium caused by pH change) and can be made after incubation at 60 °C for 48 hours and/or for 7 days.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *Geobacillus stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

T = 15-30 °C, RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place.

*ATCC® is a registered trademark of American Type Culture Collection.

Infection Control Division | **Chemical Indicators**

KPCD222-C Steam Test Pack PCD Kit

Process Challenge Device for Steam sterilization processes

Transport conditions

It is recommended to store it away from sunlight and at a temperature between 10-30 °C.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

Acceptable excursions during transport are as follows:

T = 2-10 °C and T = 30-37 °C for no longer than seven days. No restrictions for relative humidity.

Shelf life

2 years.

Packing

Presentations:

25 PCD222-C + 25 BT222 Biological Indicators per box.

Packing information: product description, storage conditions and manufacturer information.

Labeling

On product's packing: product code and description, presentation, regulation, batch number, batch number of the Self-contained Biological Indicator, manufacture and expiration date, bar code and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other relevant information

Read product's instructions for use thoroughly before use.

Precautions

Do not store the product near sterilizing agents.
Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

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