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Spare Kit for Replacement in Helix Process Challenge Devices (PCD)

RCDPCD

English

Product description

The Chemdye® Spare Kit contains a high resistance plastic tube of 2 mm internal diameter (available in several lengths). The plastic tube is arranged in a circular way (helix) tied by seals (A), it should not be released during the assembly or during the use of the device. The kit also includes a plastic seal (B), two silicone connectors (C and D) and a heat shrinkable connector (E) (Figure 1).

Assembly instructions

- 1- Introduce the A end of the plastic helix through the heat shrinkable connector (E), and then through both silicone connectors, first through the bigger one (D) and then through the smaller one (C). **IMPORTANT:** The silicone connector (C) is only necessary when replacing the plastic helix in a metal capsule device. In a plastic capsule device only use connectors (D) and (E).
- 2- Cover the A end of the plastic helix with silicone connector (C) and introduce them in the metal capsule (this connector tightens the helix to the capsule) (Figure 2). In a plastic capsule device, avoid this step and directly introduce the A end of the plastic helix in the plastic capsule.
- 3- Cover the union between the capsule and the plastic helix with the silicone connector (D).
- 4- Put the heat shrinkable connector (E) over the silicone connector (D).
- 5- Apply heat to secure the joint.
- 6- Hold the device with the provided seal (B) (Figure 2).

Additional technical information

In each cycle, the user should check that the connectors and the plastic helix are in perfect conditions. In case the plastic helix is blocked or has leakages, the complete device should be replaced.

Chemdye® Spare Kit is manufactured, calibrated and designed according to EN 285, EN 867-5 standards. Use capsules recommended by Terragene for the different types of sterilization processes.

Caution

Do not adjust the seals excessively in order to avoid choking the plastic helix.

Storage

Chemdye® Spare Kit should be stored in dry places, protected from light. Do not store close to sterilizing agents.