

Rev.0 | 11.2019



Trazanto

Start-up

Trazanto Analysis Module is incorporated to the latest version of Bionova® Traceability Software. Please download and install the latest software version before connecting the Trazanto Reader for the first time.

You can download the Bionova® Traceability Software by accessing our website www.terragene.com.ar/software.

Thank you for choosing our products!

Puesta en marcha

El Módulo de Análisis de Trazanto es incorporado en la última versión del Software de Trazabilidad Bionova®. Por favor descargue e instale la última versión de Software antes de conectar el Lector Trazanto por primera vez.

Puede descargar el Software de Trazabilidad Bionova® accediendo a nuestro sitio web www.terragene.com.ar/software.

Gracias por elegir nuestros productos!

TERRAGENE®