


Result interpretation guide

Guía de interpretación de resultados


Guia de interpretação dos resultados

Chemdye® CDSUV-1


Unprocessed | Sin procesar | Sem processar




Radiation dose | Dosis de radiación | Dose de radiação




50 mJ/cm2



100 mJ/cm2



150 mJ/cm2



200 mJ/cm2

NOTE: Screen calibration of electronic devices (PC, tablet, mobile phone), as well as environmental light, may interfere in the representation of the colors. Terragene® does not guarantee that the color you see on the screen is 100% real.

NOTA: La calibración de la pantalla de los dispositivos electrónicos (PC/tablet/teléfono móvil), como así también la luz ambiental, pueden interferir en la representación de los colores. Terragene® no garantiza que el color que ve en pantalla sea 100% real.

NOTA: A calibração da tela dos dispositivos eletrônicos (PC/tablet/celular), assim como a luz ambiente, podem interferir na representação das cores. Terragene® não garante que a cor apreciada na tela seja 100% real.